

## Handheld Pulse Oximeter

Model: SA210



- **%SpO<sub>2</sub> & Pulse Rate**  
Backlight LCD display for easy reading %SpO<sub>2</sub> & Pulse Rate
- **Personal SpO<sub>2</sub>, Pulse Rate, Memory Interval thresholds setting**
- **Shielded Design blocks Ambient Light**
- **Biocompatibility & Anti-Allergic Design**
- **Alarm Off**  
Device delivers audible warning when measurement result exceeds the SpO<sub>2</sub>/Pulse Rate thresholds previously set. User can disable the alarm function by pressing the "Alarm off" button.
- **Silent mode**  
Disables the beep sound during the measurement by pressing the "mute" button.
- **Up to 288 Memories**  
Could be programmed with 5 intervals: 0/5/10/30/60-minute ("0" means no memory setup)
- **Pulse Strength Indicator**



%SpO<sub>2</sub> & Pulse Rate



High/Low SpO<sub>2</sub>, PR & Memory Interval Setting



Shielded Design blocks Ambient Light



Biocompatibility & Anti-Allergic Design



Alarm & Sound Off



Up to 288 Memories



Backlight



Pulse Strength Indication



Protection Pouch

Optional Accessory:



Adult Probe  
PA100



Neonatal Probe  
PB100



Pediatric Probe  
PC100



Disposable Probe  
PD100



Tube Probe  
PF100

### Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring

### Specification

**Measurement Range** | SpO<sub>2</sub>: 35%~100%; PR: 30~250 bpm

**Precision** | SpO<sub>2</sub>: 70%~100%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3

**SpO<sub>2</sub> Accessory** | Adult Probe Included, Pediatric / Neonatal Probe Optional

**Power Supply** | 4 x "AA" Alkaline

Model	Qty per carton	Carton volume
SA210	10 pcs	0.039 cbm/ctn

- Storage temperature: -25°C - 70°C (-13°F - 158°F)
- Relative humidity: 15% - 90% (no condensing)
- Atmospheric Pressure: 700 to 1060 hPa