

## Fingertip Pulse Oximeter

Model: SB220



- Instant readings of SpO<sub>2</sub>, Pulse rate
- Shield Design against ambient light
- Biocompatibility protection
- Pulse Strength Indicator
- Visible & audible alarm
- Two way, two color OLED display
- String-attached



%SpO<sub>2</sub> & Pulse Rate



Shield Design against Ambient Light



Biocompatibility protection



Pulse Strength Indication



Visible/Audible Alarm



OLED Two-way Two-Color Display



String-Attached

### Specification

- Measurement range: SpO<sub>2</sub>: 35%~99%; PR: 30~250 bpm
- Precision: SpO<sub>2</sub>: 70%~99%; ±2%; 35%~69% (unspecified); Pulse rate: 30~250bpm; ±3
- Alarm: Default Value for SpO<sub>2</sub>
- Power Supply: 2"AAA" Alkaline

| Model | Qty per carton | Carton volume |
|-------|----------------|---------------|
| SB220 | 40 pcs         | 0.028 cbm/ctn |

● Storage temperature: -25°C - 70°C (-13°F - 158°F)  
Relative humidity: 15% - 90% (no condensing)  
Atmospheric Pressure: 700 to 1060 hPa

### Pulse Oximeter

- Ideal for oxygen therapy assessment
- Ideal for sleep apnea screening
- Ideal for health management
- Ideal for spot-check & continuous monitoring