

Scissors

BYDAND INSTRUMENT CATALOGUE
VERSION 11






The Supercut Advantage

All Supercut Scissors are expertly manufactured from the highest quality stainless steel. The result is a significantly sharper instrument with super fine precision capabilities, making them ideal for a variety of surgical procedures.

1. What are Supercut Scissors ?

Supercut Scissors offer unrivalled quality and cutting ability and are "sharper" than conventional scissors. The cutting blade is "razor-edged".



2. Why Supercut Scissors are better than the ordinary scissors?

A specially designed razor-sharp upper blade edge that cuts effortlessly through tissue. combined with a lower blade having micro serrations to hold and prevent tissue slippage make these scissors ideal for gripping tissue and cutting interproximally.



3. What are their other advantages?

- Supercut Scissors with razor edge + micro serrations ensure extremely controlled cutting in difficult areas during surgery. They cut their way forward more accurately.
- These scissors supersede because of their ability to cut with extreme precision right down to the tip of the blades.
- Literally speaking, the cutting edges obey the hands and mind of the surgeon most precisely.
- The cutting margins/ edges are honed to eliminate any traumatization of the tissues during cutting action, by applying special custom-made technology.

For easy identification all **Hipp** and **Armo** Supercut Scissors have black ring handles.



Contents

Braun Stadler Scissors	9	Mayo Harrington Scissors	11	Spencer Suture Scissors	15
Castroviejo Scissors	5	Mayo Scissors	10	Stevens Scissors	8
Cuticle Scissors	5	Mayo Scissors Tungsten Carbide	10	Stevens Supercut Scissors	8
Fergusson Scissors	8	Mayo Supercut Scissors	10	Strabismus Scissors	7
GoldmanN Fox Supercut Scissors	15	Metzenbaum Scissors	11	Surgical Scissors	13
Iris Scissors	6	Metzenbaum Scissors Fine	11	Surgical Scissors Tungsten Carbide	14
Iris Scissors Tungsten Carbide	6	Metzenbaum Scissors Tungsten Carbide	12	Umbilical Scissors	9
Iris Supercut Scissors	6	Metzenbaum Supercut Scissors	12	Universal Scissors	4
Kelly Scissors	7	Nail Scissors	5	Wagner Supercut Scissors	15
Kilner Scissors	7	O'Brien Scissors	8	Ward Scissors	5
Kilner Supercut Scissors	7	Reynolds Supercut Scissors	15		
Lister Bandage Scissors	4	Sims Uterine Scissors	9		



FOUR CHOICES OF QUALITY

Our selection of HIPPE, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GmbH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPE name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



LISTER BANDAGE SCISSORS

Lister Bandage Scissors are used for trimming dressings and removing bandages as well as cutting tougher material such as plaster and fabric. The blades are angled, and the lower blade is slightly longer. Available in 14.5cm, 18cm and 20cm lengths.

	HIPP	KLINI	ARMO	SAYCO
14.5cm	1.507.14	KI.507.14	A3138	280
18cm			A3139	281
20cm			A3140	



UNIVERSAL SCISSORS

Universal/Trauma scissors are used by paramedics and other emergency medical personnel to cut clothing. They have a large finger ring for easy handling and angled serrated blade for precise cutting. They are autoclavable and have coloured handles and plastic tip.

SAYCO	
19cm - LARGE	
Black	380 Black
Green	380 Green
Red	380 Red
Pink	380 Pink
16cm - SMALL	
Black	310 Black
Green	310 Green
Red	310 Red
Blue	310 Blue
Orange	310 Orange
Pink	310 Pink
Purple	310 Purple
Yellow	310 Yellow



WARD SCISSORS

Ward scissors are general purpose scissors commonly used for cutting bandages and gauze, available with or without black rubber coated handle.

SAYCO

20cm

379

20cm

379B



CUTICLE SCISSORS

Cuticle scissors are finely pointed and thin with a small blade in order to trim tiny hangnails or frayed tissue around the nail.

ARMO

STRAIGHT

9cm

A1946

CURVED

9cm

A1947



NAIL SCISSORS

Nail scissors are small scissors with a slender shaft and strong blades that are used for shaping and trimming fingernails.

ARMO

STRAIGHT


9cm

A1948

CURVED

9cm

A1949



CASTROVIEJO SCISSORS

Castroviejo scissors are double spring ophthalmic scissors commonly used during iridectomy and corneal surgery. Available in 10cm length with blunt or sharp tips.

ARMO

BLUNT/BLUNT CURVED

10cm

A3144

SHARP/SHARP CURVED

10cm

A3145



IRIS SCISSORS

Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT				
9cm	1.304.09	K1.304.09	A3150	
11cm	1.304.11	K1.304.10	A3152	299
CURVED				
9cm	1.305.09	K1.305.09	A3151	
11cm	1.305.11	K1.305.10	A3153	300



IRIS SCISSORS TUNGSTEN CARBIDE

Ward scissors are general purpose scissors commonly used for cutting bandages and gauze, available with or without black rubber coated handle.

	HIPP	KLINI	ARMO
STRAIGHT			
11cm	TC1.874.11	KTC1.874.11	A3149
CURVED			
11cm	TC1.875.11	KTC1.875.11	A3148



IRIS SUPERCUT SCISSORS

Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications.

	HIPP	KLINI	ARMO
STRAIGHT			
11cm	S1.304.11		AS3152
CURVED			
11cm	S1.305.11		AS3153



KILNER SCISSORS

Kilner scissors are used for the dissecting and cutting of delicate tissue. They are available in 13cm length with straight or curved blades. **Also available with supercut blades, please see our supercut scissors catalogue**

	HIPP	KLINI	ARMO
STRAIGHT			
13cm	1.344.13		A3155
CURVED			
13cm	1.345.13	K1.345.13	A3156



KILNER SUPERCUT SCISSORS

Kilner scissors are used for the dissecting and cutting of delicate tissue. They are available in 13cm length with straight or curved blades.

	HIPP	KLINI	ARMO
CURVED			
13cm	S1.345.13		AS3156



STRABISMUS SCISSORS

Strabismus scissors are blunt tipped scissors used for ophthalmic dissection. Available in straight or curved in 11cm length.

	HIPP	KLINI	ARMO
STRAIGHT			
11cm	1.250.11	K1.250.11	A3157
CURVED			
11cm	1.251.11	K1.251.11	A3158



KELLY SCISSORS

Kelly scissors are sharp tipped and are used to grasp and dissect fascia and soft tissue.

	HIPP	KLINI	ARMO
STRAIGHT			
16cm	1.334.16		A3159
CURVED			
16cm	1.335.16		A3160



Scissors

08

BYDAND INSTRUMENT CATALOGUE VERSION 11



STEVENS SCISSORS

Stevens Tenotomy Scissors are used for delicate dissection and cutting, usually in ophthalmic, neurological, and plastic surgery procedures. These scissors have 11cm long handles with small blades, straight or curved.

	HIPP	KLINI	ARMO
STRAIGHT			
11cm	1.270.11	K1.270.11	A3161
CURVED			
11cm	1.271.11	K1.271.11	A3162



STEVENS SUPERCUT SCISSORS

Stevens Tenotomy Scissors are used for delicate dissection and cutting, usually in ophthalmic, neurological, and plastic surgery procedures. These scissors have 11cm long handles with small blades, straight or curved.

	HIPP	KLINI	ARMO
STRAIGHT			
11.5cm	S1.270.11		AS3161
CURVED			
11.5cm	S1.271.11		AS3162



O'BRIEN SCISSORS

O'Brien Scissors are primarily used for suture removal.

	HIPP	KLINI	ARMO
9cm	1.297.09		
11.5cm			A3163



FERGUSSON SCISSORS

Ferguson Abdominal Scissors, have an angled profile and are used for multipurpose cutting and dissecting.

	HIPP	KLINI	ARMO
18cm		K1.047.18	
21cm			A3179



BRAUN STADLER SCISSORS

Braun-Stadler Episiotomy Scissors are used to perform episiotomies during labour and delivery. Available in 14.5cm and 22cm lengths with blunt tips.

	HIPP	KLINI	ARMO
14.5cm	1.217.14	KI.217.14	A3180
22cm	1.217.22		A3181



UMBILICAL SCISSORS

Umbilical scissors used for cutting the umbilical cord. These 10.5cm long scissors have wide, straight, curved blades.

	HIPP	KLINI	ARMO
10.5cm	1.200.10		A3182
12.5cm			A3183



SIMS UTERINE SCISSORS

Sims (Uterine) scissors have long blades, primarily used for deep tissue dissection. These scissors are available with blunt tips with straight or curved blades in 20cm and 23cm length.

	HIPP	KLINI	ARMO
STRAIGHT			
20cm	1.161.20		A3184
CURVED			
20cm	1.160.20		A3185
23cm			A3178



MAYO SCISSORS

Mayo Scissors have blunt ends and are generally used for cutting sutures and ligatures or trimming the skin around a wound. Mayo scissors are available in lengths from 14.5cm up to 22cm with or without tungsten carbide blades with straight or curved blades.

	HIPP	KLINI	ARMO
STRAIGHT			
14.5cm	1.030.15	K1.030.14	A3202
16cm		K1.030.15	A3204
18cm	1.030.17	K1.030.17	A3207
20cm	1.030.19	K1.034.19	A3208
22cm	1.030.23	K1.030.23	A3215
CURVED			
14.5cm	1.031.15	K1.031.14	A3203
16cm		K1.031.15	A3205
18cm	1.031.17	K1.031.17	A3209
20cm	1.031.19		A3211
22cm	1.031.23	K1.031.23	A3214



MAYO SCISSORS TUNGSTEN CARBIDE

Mayo Scissors have blunt ends and are generally used for cutting sutures and ligatures or trimming the skin around a wound. Mayo scissors are available in lengths from 14.5cm up to 22cm with or without tungsten carbide blades with straight or curved blades.

	HIPP	KLINI	ARMO
STRAIGHT			
14.5cm	TC1.810.14	KTC1.810.14	A3200
18cm	TC1.810.17	KTC1.810.17	A3206
22cm	TC1.810.23	KTC1.810.23	A3212
CURVED			
14.5cm	TC1.811.14	KTC1.811.14	A3201
17cm	TC1.811.17	KTC1.811.17	A3210
22cm	TC1.811.23	KTC1.811.23	A3213



MAYO SUPERCUT SCISSORS

Mayo Scissors have blunt ends and are generally used for cutting sutures and ligatures or trimming the skin around a wound.

	HIPP	KLINI	ARMO
STRAIGHT			
14cm	S1.030.14		A53202
17cm	S1.030.17		A53207
20cm	S1.030.20		A53208
CURVED			
14cm	S1.031.14		A53203
17cm	S1.031.17		A53209



MAYO HARRINGTON SCISSORS

Mayo Harrington Dissecting Scissors have flat sharp blades with blunt tips and are available in straight or curved in a 22cm length.

HIPP

STRAIGHT	
22cm	1.030.22
CURVED	
22cm	1.031.22



METZENBAUM SCISSORS

Metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. These scissors feature a more delicate profile than the standard metzenbaum scissor.

HIPP

KLINI

ARMO

STRAIGHT			
12cm			A3249
14cm	1.050.14	K1.050.14	A3253
16cm	1.070.15		A3255
18cm	1.050.18	K1.050.18	A3257
20cm	1.050.20	K1.050.20	A3262
22cm	1.050.23		A3263
CURVED			
12cm	1.051.12		A3250
14cm	1.051.14	K1.051.14	A3254
16cm	1.071.15		A3256
18cm	1.051.18	K1.051.18	A3258
20cm	1.051.20	K1.051.20	A3264
22cm	1.051.23		A3265
28cm			A3248



METZENBAUM SCISSORS FINE


Metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than Mayo scissors.

HIPP

KLINI

ARMO



FINE TIP STRAIGHT	
15cm	A3246
FINE TIP CURVED	
15cm	A3247



METZENBAUM SCISSORS TUNGSTEN CARBIDE

Tungsten carbide metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than mayo scissors. Having the blades infused with tungsten carbide allows for greater durability of the scissor.

	HIPP	KLINI	ARMO
STRAIGHT			
14cm	TC1.830.14	KTC1.830.14	A3251
18cm	TC1.830.18	KTC1.830.18	A3266
20cm	TC1.830.20	KTC1.830.20	A3267
22cm	TC1.830.23		A3260
CURVED			
14cm	TC1.831.14	KTC1.831.14	A3252
18cm	TC1.831.18	KTC1.831.18	A3259
20cm	TC1.831.20	KTC1.831.20	A3268
22cm	TC1.831.23		A3261



sharp/blunt

blunt/blunt

METZENBAUM SUPERCUT SCISSORS

Supercut metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than Mayo scissors.

STRAIGHT			
	HIPP	KLINI	ARMO
SHARP/BLUNT			
14cm	S1.052.14		AS3166
18cm	S1.062.18		AS3167
BLUNT/BLUNT			
11cm	S1.050.11		
14cm	S1.050.14		AS3253
18cm	S1.060.18		AS3257
CURVED			
	HIPP	KLINI	ARMO
SHARP/BLUNT			
14cm	S1.053.14		AS3168
18cm	S1.063.18		AS3169
BLUNT/BLUNT			
11cm	S1.051.11		AS3250
14cm	S1.051.14		AS3254
18cm	S1.061.18		AS3258



SURGICAL SCISSORS

Surgical or operating scissors are used for general cutting and dissecting of tissue, bandages and gauze. They are available with blunt or sharp blades or a combination of both. Surgical scissors are available in lengths from 11cm up to 20cm with straight or curved blades.




sharp/sharp



sharp/blunt




blunt/blunt



sharp/probe



sharp/sharp



sharp/blunt



blunt/blunt

STRAIGHT

	HIPP	KLINI	ARMO
SHARP/SHARP			
11cm		K1.004.11	A3296
13cm	1.004.13	K1.004.13	A3302
14cm		K1.004.14	A3290
16cm	1.004.16	K1.004.16	A3320
18cm	1.004.18		A3323
20cm			A3346
SHARP/BLUNT			
11cm		K1.002.11	A3297
13cm	1.002.13	K1.002.13	A3301
14cm	1.002.14	K1.002.14	A3292
16cm	1.002.16	K1.002.16	A3319
18cm	1.002.18		A3339
BLUNT/BLUNT			
11cm			A3298
13cm	1.000.13	K1.000.13	A3300
14cm			A3294
16cm	1.000.16	K1.000.16	A3318
18cm	1.000.18		
SHARP/PROBE			
13cm			A3306

CURVED

	HIPP	KLINI	ARMO
SHARP/SHARP			
11cm		K1.005.11	A3295
13cm	1.005.13	K1.005.13	A3305
14cm		K1.005.14	A3289
16cm	1.005.16	K1.005.16	A3317
18cm	1.005.18		A3332
20cm			
SHARP/BLUNT			
11cm			
13cm	1.003.13	K1.003.13	A3304
14cm		K1.003.14	A3291
16cm	1.003.16	K1.003.16	A3316
18cm	1.003.18		A3331
BLUNT/BLUNT			
11cm			A3299
13cm	1.001.13	K1.001.13	A3303
14cm			A3293
16cm	1.001.16	K1.001.16	A3315
18cm	1.001.18		A3330



SURGICAL SCISSORS TUNGSTEN CARBIDE

Tungsten carbide surgical or operating scissors are used for general cutting and dissecting of tissue, bandages and gauze. They are available with blunt or sharp blades or a combination of both. Having the blades infused with tungsten carbide allows for greater durability of the scissor.

	STRAIGHT		
	HIPP	KLINI	ARMO
SHARP/SHARP			
14cm	TC1.804.14	KTC1.804.14	A3307
16cm	TC1.804.16		A3324
18cm	TC1.804.18		A3333
SHARP/BLUNT			
14cm	TC1.802.14	KTC1.802.14	A3308
16cm	TC1.802.16		A3325
18cm			A3335
BLUNT/BLUNT			
14cm	TC1.800.14	KTC1.800.14	A3309
16cm	TC1.800.16		A3326
18cm			A3337

	CURVED		
	HIPP	KLINI	ARMO
SHARP/SHARP			
14cm	TC1.805.14		A3310
16cm	TC1.805.16		A3327
18cm			A3334
SHARP/BLUNT			
14cm	TC1.803.14		A3311
16cm	TC1.803.16		A3328
18cm			A3336
BLUNT/BLUNT			
14cm	TC1.801.14		A3312
16cm	TC1.801.16		A3329
18cm			A3338

SPENCER SUTURE SCISSORS

Spencer Scissors are primarily used for suture removal. These scissors have a hook shaped tip on one blade that slides under sutures to slightly lift them before cutting. This hook also prevents the suture from slipping off the blade prior to cutting.

	HIPP	KLINI	ARMO
9cm	1.372.09	K1.372.09	A3400
12cm	1.374.12	K1.372.12	A3401
14cm			A3402



GOLDMANN FOX SUPERCUT SCISSORS

Goldmann Fox Scissors are used in dental procedures to trim gingival tissue. They are also used in plastic surgery or hemorrhoidectomy. They have curved, bevelled blades, one of which is serrated, that taper into sharp, fine tips. Available in 13cm length with Supercut blades

	HIPP	KLINI	ARMO
CURVED			
13cm	S1.487.13Z		AS3165



REYNOLDS SUPERCUT SCISSORS

Reynolds scissors are characterised by having narrow, tenotomy style tips, generally used for cutting tissue as well as dissection. Available in 15cm length with curved supercut blades.

	HIPP	KLINI	ARMO
CURVED			
15cm	S1.325.15		AS3147



WAGNER SUPERCUT SCISSORS

Wagner Scissors are used in procedures requiring delicate or precise cutting, particularly plastic surgery. Available in 12cm length with straight or curved Supercut blades.

	HIPP	KLINI	ARMO
STRAIGHT			
12cm	S1.304.12		AS3163
CURVED			
12cm	S1.305.12		AS3164



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au