

Forceps

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Contents (by kind)	2
Contents (alphabetically ordered)	2

Aural Dressing Forceps	12
Catheter Forceps	13
Dental Forceps	14
Dissecting Forceps	15
Dressing Forceps	4
Entomological/Jewellers Forceps	15
Eye Forceps	16
First Aid & Splinter Forceps	10
Intestinal Clamps/Forceps	20
Sponge Holding Forceps	21
Steriliser Forceps	22
Tissue Forceps	6
Uterine Forceps	22
Vasectomy Forceps	24



Contents

Adson	6	Entomological/Jewellers Forceps	15	Perry Dental, Serrated Tips	14
Adson Brown	7	Feilchenfeld	10	Potts Smith	9
Adson-tungsten carbide	7	First Aid/splinter	10	Potts Smith, Serrated Tips	5
Adson, Serrated Tips	4	Flexible Alligator	13	Rampley	21
Alderkreutz	6	Foerster	21	Roberts	19
Allis	8	Forceps with Magnifier	10	Schroeder-Braun, Tenaculum	23
Ayer Chalazion	16	Gillies	8	Schroeder, Tenaculum	22
Babcock	8	Graefe Iris	18	Semken	9
Barraquer Cilia	17	Hartman Polypus, Serrated Jaws	12	Semken, Serrated Tips5	
Beer Cilia	16	Hartman Standard	12	Skene	22
Bishop Harman	19	Hunter	11	Standard	6
Block End	4	Jansen Aural Dressing, Serrated Ends	13	Stiegiltz	11
Bonney	8	Lambert Chalazion	16	Teale, Vulsellum	23
Bozeman, Dressing	23	Lane	9	Tilley Aural Dressing, Serrated Ends	12
Castroviejo Corneal	19	Lester, 1x2 Teeth	18	Treves	7
Castroviejo, 1x2 Teeth	18	Lister Sinus	21	Troeltsch Aural Dressing, Serrated Ends	13
Cheatle	22	London College, Serrated Tips	14	Tubing Clamp	20
Cilia, Slant Edge	17	Magill	13	Turnover End	4
Colibri	19	Maier Polypus	21	Vasectomy	24
Desmarres Chalazion	16	Mayo-Robson/ Doyen	20	Vasectomy Towel Clamps	24
Dissecting, Pointed	15	McGee Micro Ear, Micro Serrated Jaws	12		
Douglas Cilia	17	McIndoe, Serrated Tips	5		
Doyen Intestinal Clamp	20	Pennington	20		

Contents (by kind)

Aural Dressing Forceps		Barraquer Cilia	17	Rampley	21
Flexible Alligator	13	Beer Cilia	16	Steriliser Forceps	
Hartman Polypus, Serrated Jaws	12	Bishop Harman	19	Cheatle	22
Hartman Standard	12	Castroviejo, 1X2 Teeth	18	Tissue Forceps	
Jansen Aural Dressing, Serrated Ends	13	Castroviejo Corneal	19	Adson	6
McGee Micro Ear, Micro Serrated Jaws	12	Cilia, Slant Edge	17	Adson Brown	7
Tilley Aural Dressing, Serrated Ends	12	Colibri	19	Adson-tungsten carbide	7
Troeltsch Aural Dressing, Serrated Ends	13	Desmarres Chalazion	16	Alderkreutz	6
Catheter Forceps		Douglas Cilia	17	Allis	8
Magill	13	Graefe Iris	18	Babcock	8
Dental Forceps		Lambert Chalazion	16	Bonney	8
London College, Serrated Tips	14	Lester, 1X2 Teeth	18	Gillies	8
Perry Dental, Serrated Tips	14	Roberts	19	Lane	9
Dissecting Forceps		First Aid & Splinter Forceps		Potts Smith	9
Dissecting, Pointed	15	Feilchenfeld	10	Semken	9
Dressing Forceps		First Aid/splinter	10	Standard	6
Adson, Serrated Tips	4	Forceps with Magnifier	10	Treves	7
Block End	4	Hunter	11	Uterine Forceps	
McIndoe, Serrated Tips	5	Stiegiltz	11	Bozeman, Dressing	23
Potts Smith, Serrated Tips	5	Intestinal Clamps/Forceps		Schroeder-Braun, Tenaculum	23
Semken, Serrated Tips	5	Doyen Intestinal Clamp	20	Schroeder, Tenaculum	22
Turnover End	4	Mayo-Robson/ Doyen	20	Skene	22
Entomological/Jewellers Forceps		Pennington	20	Teale, Vulsellum	23
Jewellers	15	Tubing Clamp	20	Vasectomy Forceps	
Eye Forceps		Sponge Holding Forceps		Vasectomy	24
Ayer Chalazion	16	Foerster	21	Vasectomy Towel Clamps	24
		Lister Sinus	21		
		Maier Polypus	21		



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

 <p>ANTON HIPPI Instruments & Implants</p>		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
 <p>KLINI™ instruments by anton hipp GERMANY</p>		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
 <p>ARMO surgical instruments</p>		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
 <p>Sayco first aid instruments</p>		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



Dressing Forceps

BLOCK END

Dressing Forceps are designed to hold gauze and other dressings. They are for general use including grasping dressings, suture material, tissue or foreign bodies in wounds.

	HIPP	KLINI	ARMO	SAYCO
BLUNT				
11cm			A2404	
13cm	3.000.13	K3.000.13	A2402	109
14cm			A2138	
16cm	3.000.16	K3.000.16	A2139	110
18cm		K3.000.18	A2409	111
20cm			A2140	
25cm			A2141	
DELICATE				
13cm		K3.004.13	A2403	

TURNOVER END

Dressing Forceps are designed to hold gauze and other dressings. They are for general use including grasping dressings, suture material, tissue or foreign bodies in wounds.

	HIPP	KLINI	ARMO	SAYCO
BLUNT				
13cm				103
15cm				104

ADSON, SERRATED TIPS

Adson Dressing Forceps are used for holding materials such as cotton and gauze, changing dressings, or packing wounds. They have a wide thumb grasp for improved control.

	HIPP	KLINI	ARMO
STANDARD POINT			
12cm	3.040.12	K3.040.12	A2401
15cm			A2407
MICRO POINT			
12cm	3.044.12	K3.044.12	A2400



Dressing Forceps

MCINDOE, SERRATED TIPS

McIndoe Forceps are used for holding materials such as cotton and gauze, changing dressings, or packing wounds. These forceps have a stop peg to ensure that when used for holding tissue the grip is not crushing.

	HIPP	KLINI	ARMO	SAYCO
15cm	3.050.15	K3.050.15	A2415	101

POTTS SMITH, SERRATED TIPS

Potts-Smith Dressing Forceps are used for applying and changing of dressings. They may also be used for dissection of tissue.


	HIPP	KLINI	ARMO
STRAIGHT			
18cm	3.080.18	K3.080.18	A2608
21cm			A2136
CURVED			
18cm	3.081.18		

SEMKEN, SERRATED TIPS

Semken Dressing Forceps are most commonly used for applying and changing of dressings. They have narrow tips that may be beneficial for neurosurgery or other procedures performed in hard to reach places. They may also be used for dissection of fine or single-layer tissue. These forceps have a stop peg to ensure that when used for holding tissue the grip is not crushing.

	HIPP	KLINI	ARMO
12.5cm	3.054.12	K3.054.12	A2610
15cm	3.054.15		A2611





Tissue Forceps

STANDARD

Tissue Forceps are used for grasping, holding, and manipulating various types of tissue. Thumb tissue forceps are toothed to grip the tissue more firmly than the equivalent serrated tip forceps.

	HIPP	KLINI	ARMO
1X2 TEETH			
11cm			A2606
13cm			A2615
14cm	3.100.14	K3.100.14	A2616
16cm	3.100.16	K3.100.16	A2617
18cm			A2624
20cm	3.100.20		A2614
25cm	3.100.25		A2618
2X3 TEETH			
12.5cm	3.102.13	K3.102.13	A2621
14cm	3.102.14		A2622
16cm	3.102.16		A2623

ALDERKREUTZ

Alderkreutz tissue forceps are commonly used for neurosurgery and cardiovascular surgery due to the wider tips.

	HIPP	KLINI	ARMO
4X5 TEETH			
15cm	3.238.15		




ADSON

Adson Tissue Forceps are used for holding and manipulating delicate tissue.

	HIPP	KLINI	ARMO
1X2 STANDARD			
12cm	3.150.12	K3.150.12	A2600
15cm			A2592
1X2 TYING PLATFORM			
12cm	3.154.12		
1X2 MICRO			
12cm	3.158.12	K3.158.12	A2559
15cm			A2591



Tissue Forceps



ADSON-TUNGSTEN CARBIDE

Adson Tissue Forceps are used for holding and manipulating delicate tissue.

	HIPP	KLINI	ARMO
1X2 STANDARD			
12cm			A2601

ADSON BROWN

Adson Brown Tissue Forceps are used for holding and manipulating delicate tissue. They have jaws with 9x9 serrated teeth running horizontally down the tips of the forcep.

	HIPP	KLINI	ARMO
12cm	3.160.12	K3.160.12	A2554

TREVES

Treves tissue forceps are characterised by having a tip with interlocking teeth making them suitable for holding skin or dense tissue.

	HIPP	KLINI	ARMO	SAYCO
1X2 TEETH				
13cm			A2630	147
15cm				151



ALLIS

Allis Tissue Forceps are ratcheted forceps used for grasping organs and slippery or dense tissue during surgery.


	HIPP	KLINI	ARMO
4X5 TEETH			
15cm		K2.392.15	A2605
18cm			A2603
14cm			A2602



BABCOCK

Babcock Tissue Forceps are ratcheted, non-perforating forceps used to grasp delicate tissue.

	HIPP	KLINI	ARMO
16cm			A2620
14cm			A2632



BONNEY

Bonney tissue forceps are used for grasping heavy tissue or bone.

	HIPP	KLINI	ARMO
1X2 TEETH			
18cm		K3.222.18	A2619






GILLIES

Gillies Tissue Forceps are used in small surgical sites to grasp tissue. They have fine, narrow jaws that are commonly used in plastic surgery.

	HIPP	KLINI	ARMO	SAYCO
1X2 TEETH				
15cm	3.166.15	K3.166.15	A2625	150



Tissue Forceps



LANE

Lane Tissue Forceps are used for holding tough tissue. The handle is smooth and features a slight thumb/finger indentation for greater control.

	HIPP	KLINI	ARMO
1X2 TEETH			
14.5cm		K3.246.14	A2607

POTTS SMITH

Potts Smith tissue forceps are used for grasping fine or delicate tissue. They have a stop peg to prevent an overly harsh grasp that may crush the tissue.

	HIPP	KLINI	ARMO
1X2 TEETH			
18cm	3.200.18	K3.200.18	A2609

SEMKEN

Semken Tissue Forceps are lightweight, narrow forceps used for gripping fine membranes or single layer tissues.

	HIPP	KLINI	ARMO
1X2 TEETH			
12.5cm	3.170.12	K3.170.12	A2612
15cm	3.170.15	K3.170.15	A2613



First Aid & Splinter Forceps

FEILCHENFELD

Feichenfeldt splinter forceps are multi-purpose with sharp tips used for removing small foreign objects. Available in 9cm and 11.5cm lengths.

	HIPP	KLINI	ARMO	SAYCO
9cm	3.370.09	K3.370.09	A2550	
11.5cm	3.370.11	K3.370.11	A2551	138



FIRST AID/SPLINTER

First Aid/Splinter forceps can be recognized by their serrated handles that taper to fine tips. These forceps are suitable for removing splinters or other minor medical uses. Available in 9cm length with pointed, flat or slanted tips and 12.5cm with pointed tips.

	HIPP	KLINI	ARMO	SAYCO
POINTED				
9cm			A2519	
12.5cm			A2520	F830
FLAT				
9cm				600
SLANTED				
9cm				599



FORCEPS WITH MAGNIFIER

These forceps with magnifier lens aids the precise location of the splinter. Available in 10cm length straight or curved.

	HIPP	KLINI	ARMO	SAYCO
FLAT				
10cm				206
CURVED				
10cm				208




First Aid & Splinter Forceps

HUNTER

Hunter Forceps are characterised by having an arrowhead pointed serrated tip and are used for the removal of splinters and other foreign bodies or to hold and manoeuvre dressings or soft tissue during procedures.


	HIPP	KLINI	ARMO	SAYCO
STRAIGHT NO PIN				
10cm	3.378.10	K3.378.10		
STRAIGHT WITH PIN				
10cm	3.374.10	K3.374.10	A2552	
CURVED WITH PIN				
10cm	3.375.10		A2555	




STIEGILTZ

Stieglitz Forceps are ratcheted forceps with an arrowhead pointed serrated tip and are usually used in dental procedures.

	HIPP	KLINI	ARMO	SAYCO
15cm	60.058.13	K3.384.13	A2553	





Aural Dressing Forceps

HARTMAN STANDARD

The Hartmann alligator forceps or Hartmann foreign body forceps are medical forceps for removing foreign bodies. It is used in addition to surgery mainly in ENT.

	HIPP	KLINI	ARMO	SAYCO
14cm			A2209	



HARTMAN POLYPUS, SERRATED JAWS

Hartman Polypus Forceps are used for ENT applications where arteries or objects in the aural cavity need to be clamped, grasped or removed. They have serrated jaws.

	HIPP	KLINI	ARMO	SAYCO
8cm			A2200	
11cm			A2201	
14cm			A2202	
1X2 TEETH				
8cm			A2205	



MCGEE MICRO EAR, MICRO SERRATED JAWS

McGee Micro Ear Forceps are right-angled forceps with serrated jaws used in procedures to finely excise tissues in the narrow spaces of the aural cavity.

	HIPP	KLINI	ARMO	SAYCO
8cm			A2206	
11cm			A2207	




TILLEY AURAL DRESSING, SERRATED ENDS

Tilley Aural Dressing Forceps have serrated jaws and are designed for the removal or insertion of nasal dressings or packs.

	HIPP	KLINI	ARMO	SAYCO
14cm	26.601.13		A2204	




Aural Dressing Forceps



FLEXIBLE ALLIGATOR

The Flexible Alligator forcep has a flexible shaft that allows it to be fed through Veterinary Oscopes.


	HIPP	KLINI	ARMO	SAYCO
14cm			A2215	



TROELTSCH AURAL DRESSING, SERRATED ENDS

Troeltsch Forceps are commonly used in ENT procedures. These distinctive serrated forceps are bent at an 80-degree angle and are particularly useful for packing and dressing in ENT procedures.


	HIPP	KLINI	ARMO	SAYCO
13cm	3.429.13	K3.429.13	A2198	
18cm			A2185	




JANSEN AURAL DRESSING, SERRATED ENDS

Jansen dressing forceps are commonly used in nasal procedures. The bayonet style handle and long serrated jaws assist in an unobstructed view when working in the nasal cavity.


	HIPP	KLINI	ARMO	SAYCO
20cm			A2203	



Infant
1/1



Child
1/1



Adult
1/1

Catheter Forceps

MAGILL

Magill Catheter Forceps have fenestrated tips and are used for grasping endotracheal tubes during intubation and removing foreign bodies from the airway. They can also be used to grasp other types of catheters during insertion.

	HIPP	KLINI	ARMO	SAYCO
INFANT				
15cm			A2192	
CHILD				
20cm			A2193	
ADULT				
25cm			A2194	



Dental Forceps

PERRY DENTAL, SERRATED TIPS

Perry Dental Forceps are curved serrated forceps with a location peg used for cotton wad placement and removal.

	HIPP	KLINI	ARMO	SAYCO
CURVED				
13cm			A2295	



LONDON COLLEGE, SERRATED TIPS

London College Forceps are used to grasp or transfer materials in and out of the oral cavity. They have serrated tips and a locking device.

	HIPP	KLINI	ARMO	SAYCO
ANGLED WITH LOCK				
15cm			A2296	



Dissecting Forceps



DISSECTING, POINTED

Dissecting Forceps are used to handle tissue and other materials and to manipulate needles and other instruments whilst operating.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT				
12.5cm			A2370	095
CURVED				
13cm				462

Entomological/Jewellers Forceps



JEWELLERS

Jewellers Forceps are used in microsurgical procedures. The micro-fine tips that can handle extremely fine tissues and materials that standard instruments may damage.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT NO.3 SUPER FINE				
12.5cm			A2524	
STRAIGHT NO.5 SUPER FINE				
11cm			A2525	
CURVED NO.7 SUPER FINE				
13cm			A2527	462



AYER CHALAZION

Ayer Chalazion Forceps are ophthalmic forceps used to invert the eyelid.

	HIPP	KLINI	ARMO	SAYCO
9cm			A2189	



LAMBERT CHALAZION

Lambert Chalazion Forceps are ophthalmic forceps used to invert the eyelid.

	HIPP	KLINI	ARMO	SAYCO
9cm	3.796.09		A2190	



DESMARRES CHALAZION

Desmarres Chalazion Forceps are ophthalmic forceps used to invert the eyelid.

	HIPP	KLINI	ARMO	SAYCO
NO.1				
9cm			A2516	
NO.2				
9cm			A2517	



No.1
1/1


No.2
1/1



BEER CILIA

Beer Cilia Forceps are used for grasping and holding cilia/eyelashes. The tips are rounded with a raised centre.

	HIPP	KLINI	ARMO	SAYCO
9cm	3.414.09		A2191	





Eye Forceps

DOUGLAS CILIA

Cilia Forceps are used for removing cilia/eyelashes, foreign bodies and sutures around the eye and eyelid.

	HIPP	KLINI	ARMO	SAYCO
9cm			A2188	

CILIA, SLANT EDGE

Cilia Forceps are used for removing cilia/eyelashes, foreign bodies and sutures around the eye and eyelid.

	HIPP	KLINI	ARMO	SAYCO
9cm			A2186	

BARRAQUER CILIA

Barraquer Cilia Eye Forceps are an ophthalmologic instrument used for grasping or manipulating eyelashes.

	HIPP	KLINI	ARMO	SAYCO
9cm			A2187	







Figure 1




Figure 2



Figure 3




Figure 4







Figure 5



Eye Forceps

GRAEFE IRIS

Graefe Iris Forceps are used during ophthalmic procedures to grasp and hold tissue.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT - NO PIN (figure 1)				
7cm			A2508	
STRAIGHT - WITH PIN (figure 1)				
10cm	3.644.10	K3.300.10	A2511	
STRAIGHT - WITH PIN, 1X2 TEETH (figure 4)				
10cm	K3.664.10	A2507		
HALF CURVED - NO PIN (figure 2)				
7cm			A2509	
HALF CURVED - WITH PIN (figure 2)				
10cm	K3.301.10	A2512		
HALF CURVED - WITH PIN, 1X2 TEETH (figure 5)				
10cm	K3.669.10	A2505		
FULL CURVED - WITH PIN (figure 3)				
10cm	3.655.10		A2506	

LESTER, 1X2 TEETH

Lester Forceps are used during ophthalmic procedures to grasp and hold tissue.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT - WITH PIN				
9cm			A2510	
CURVED - NO PIN				
7cm			A2513	

CASTROVIEJO, 1X2 TEETH

Castroviejo Forceps are used in ophthalmic and plastic surgery procedures near the eyes. The forceps are used to grasp tissue and sutures.

	HIPP	KLINI	ARMO	SAYCO
9cm			A2514	



Eye Forceps

CASTROVIEJO CORNEAL FORCEPS

Castroviejo Corneal forceps are extremely delicate tissue forceps with a tying platform and a micro jaw.

	HIPP	KLINI	ARMO	SAYCO
10.5cm			A2182	

BISHOP HARMAN

The Bishop Harman Forceps are used during surgeries to grasp tissues in skin grafting or removing affected tissues from the body, they are designed to be used in various ophthalmologic procedures.

	HIPP	KLINI	ARMO	SAYCO
8cm			A2502	

ROBERTS

The Roberts Forcep is used in conjunction with a needleholder to tie knots in fine materials.

	HIPP	KLINI	ARMO	SAYCO
10cm			A2184	

COLIBRI

The Colibri forcep is designed to perform several functions. The tooth at the tip is ideal for holding the corneal or scleral edge when passing a suture and it can also be used for grasping iris tissue.

	HIPP	KLINI	ARMO	SAYCO
7.5cm			A2183	





Intestinal Clamps/ Forceps

PENNINGTON

Pennington Forceps are used for holding tissue that needs to be sutured or for haemostasis.

	HIPP	KLINI	ARMO	SAYCO
15cm			A2535	



TUBING CLAMP

Tubing Forceps are ratcheted forceps are used to control moderate to heavy flow of fluids or blood. They are a heavyweight instrument and can be used with heavy tubing or vessels.

	HIPP	KLINI	ARMO	SAYCO
15cm			A2530	



DOYEN INTESTINAL CLAMP

The Doyen Intestinal clamps have longitudinal serrations and a non crushing clamp, are a lightweight instrument and are typically used to grasp large organs. This instrument is available in both straight and curved and the serrations on the jaw run diagonally across the jaws.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT				
16.5cm		K33.108.16		
21cm		K33.114.20		
23cm		K33.115.23		
CURVED				
16.5cm		K33.109.16		
21cm		K33.115.20		
23cm		K33.115.23		



MAYO-ROBSON/ DOYEN

The Mayo-Robson/Doyen Intestinal forceps are used to gently but firmly grip onto tissue. It is commonly used for anastomosis or holding the bowel during intestinal procedures. The serrations on the jaw of the forceps run parallel along the length of the jaws.

	HIPP	KLINI	ARMO	SAYCO
STRAIGHT				
17cm				A2720
21cm				A2722
24cm				A2724
CURVED				
17cm				A2721
21cm				A2723
24cm				A2725



Sponge Holding Forceps



FOERSTER

Foerster Sponge Forceps are used in surgical procedures to hold gauze squares that will absorb excess fluids and blood from the surgical site. The rounded, looped tips are serrated.


	HIPP	KLINI	ARMO	SAYCO
18cm			A2538	
20cm			A2542	
25cm			A2537	



RAMPLEY

Rampley Sponge Forceps are used to hold sponges and swabs. The jaws are elongated and contain serrated tips to hold swabs or sponges.



	HIPP	KLINI	ARMO	SAYCO
18cm	2.318.18		A2539	
24cm	2.318.25	K2.318.25	A2540	



LISTER SINUS

Lister Sinus Forceps are used for inserting or removing packing in the sinus cavity.

	HIPP	KLINI	ARMO	SAYCO
13cm			A2521	
15cm			A2522	



MAIER POLYPUS

The Maier Polypus Sponge Holding Forcep is used for the placement, holding and removal of dressings. It's serrated tips allow a secure grips on swabs, gauze and other dressing materials.

	HIPP	KLINI	ARMO	SAYCO
25cm	2.286.25			



Steriliser Forceps

CHEATLE

Cheatle Sterilizer Forceps are used to remove sterilised instruments from the sterilisation container.

	HIPP	KLINI	ARMO	SAYCO
LIGHT PATTERN				
27cm	2.331.27		A2541	



Uterine Forceps

SKENE

Skene Forceps are ratcheted curved forceps and are used to grasp and steady the cervix and the uterus during gynaecological surgery. The tips have hooked single toothed jaws.

	HIPP	KLINI	ARMO	SAYCO
24cm	2.485.24	K2.485.24	A2701	



SCHROEDER, TENACULUM

Schroeder Tenaculum Forceps are ratcheted forceps with use in obstetric and gynaecological procedures.

	HIPP	KLINI	ARMO	SAYCO
2X2 TEETH				
25cm			A2702	

Uterine Forceps

SCHROEDER-BRAUN, TENACULUM

Schroeder Tenaculum Forceps are ratcheted forceps with use in obstetric and gynaecological procedures.

	HIPP	KLINI	ARMO	SAYCO
25cm	2.480.25	K2.480.25	A2700	

BOZEMAN, DRESSING


Bozeman Uterine Dressing Forceps are used in gynaecological procedures. The shanks are slightly curved to provide greater visibility when used in the uterine cavity.

	HIPP	KLINI	ARMO	SAYCO
SERRATED, CURVED				
26cm	2.301.26	K2.301.26	A2628	

TEALE, VULSELLUM

Teale Vulsellum Forceps are curved ratcheted forceps used in obstetric and gynaecological surgery.

	HIPP	KLINI	ARMO	SAYCO
CURVED				
25cm		K2.481.23	A2626	





Vasectomy Forceps

VASECTOMY TOWEL CLAMPS

Vasectomy Towel Clamps are ratcheted forceps used in non-suture vasectomy procedures.

	HIPP	KLINI	ARMO	SAYCO
3MM TIP				
14cm			A5171	



VASECTOMY

Vasectomy Forceps are ratcheted forceps with sharp tips used to pierce the scrotal skin in vasectomy procedures.

	HIPP	KLINI	ARMO	SAYCO
CURVED				
12cm			A5172	





19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au