

Eye & Nasal Specula

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Eye Specula

Alfonso	4
Bangerter	5
Barraquer Wire	3
Barraquer Wire (Solid)	3
Castroviejo	6
Cook	5
Knapp-Clark	6
Kratz-Barraquer Heavy (1mm) Wire	4
Kratz-Barraquer Wire	3
Lester-Burch	5
Lieberman V-Wire	6
Mellinger	5
Murdoch	4
Sauer	4
Weiss	6
Wire Specula - Infant	3

Nasal Specula

Thudichum	7
-----------	---




FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>




Eye Specula



BARRAQUER WIRE

Barraquer Wire Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery.


	HIPP	KLINI	ARMO
ADULT			
15mm blades	100 4014		



BARRAQUER WIRE (SOLID)

Barraquer Wire Solid Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery.


	HIPP	KLINI	ARMO
ADULT			
15mm blades	106 4015		



WIRE SPECULA - INFANT

Wire Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery.


	HIPP	KLINI	ARMO
INFANT			
5mm blades	108 3604		



KRATZ-BARRAQUER WIRE

Kratz-Barraquer Wire Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open usually during cataract removal surgery.

	HIPP	KLINI	ARMO
ADULT			
15mm blades	109 4015		




Eye Specula

KRATZ-BARRAQUER HEAVY (1MM) WIRE

Kratz-Barraquer Wire Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open usually during cataract removal surgery. Heavy duty wire provides greater force than standard wire.


	HIPP	KLINI	ARMO
ADULT			
15mm blades	110 4015		



ALFONSO

Alfonso Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during paediatric surgery with the assistance of thumb rests on the frame.

	HIPP	KLINI	ARMO
8 x 5mm blades	128 4608		



MURDOCH

Murdoch eye specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved retractor blades which are mounted back to back on a bar along which one of them can be moved to part the eyelids.

	HIPP	KLINI	ARMO
ADULT, LEFT			
15mm blades	138 5010		
ADULT, RIGHT			
15mm blades	138 5001		
CHILD, LEFT			
10mm blades	138 4410		
CHILD, RIGHT			
10mm blades	138 4401		



SAUER

Sauer Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during paediatric surgery with the assistance of sprung hinged blades.

	HIPP	KLINI	ARMO
11mm blades	162 3010		



Eye Specula



BANGERTER

Bangerter Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved retractor blades which are mounted back to back on a bar along which one of them can be moved to part the eyelids. A locking nut is fitted as well.

	HIPP	KLINI	ARMO
ADULT			
15mm blades	140 5515		

COOK

Cook Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved retractor blades which are mounted back to back on a bar along which one of them can be moved to part the eyelids. A locking nut is fitted as well.

	HIPP	KLINI	ARMO
CHILD			
11mm blades	136 4011		



MELLINGER

Mellinger Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved retractor blades which are mounted back to back on a spring-loaded bar along which one of them can be moved to part the eyelids.


	HIPP	KLINI	ARMO
CHILD			
10mm	134 5011		



LESTER-BURCH

Lester-Burch Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved retractor blades which are mounted on sprung pivoting arms with thumb grips and a locking nut is fitted as well.

	HIPP	KLINI	ARMO
ADULT			
15mm blades	165 8522		



Eye Specula

LIEBERMAN V-WIRE

Lieberman V-wire Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved open wire blades which are mounted with stepless adjustable spreading arms.

	HIPP	KLINI	ARMO
ADULT			
15mm blades	164 6215V		

CASTROVIEJO

Castroviejo Eye Specula is an ophthalmic instrument used to retract and hold the eyelid open during surgery. It is characterised by having curved open wire blades which are mounted with stepless adjustable spreading arms.

	HIPP	KLINI	ARMO
ADULT			
15mm blades			A5300

WEISS

The Weiss eye speculum is an ophthalmic instrument used to keep the eyelids apart during surgical intervention such as glaucoma surgery, corneal transplant, laser eye surgery or a general ocular examination.

	HIPP	KLINI	ARMO
ADULT			
	142 7010		
CHILD			
	142 9020		

KNAPP-CLARK

The Knapp-Clark eye speculum is an ophthalmic instrument used to keep the eyelids apart during surgical intervention such as glaucoma surgery, corneal transplant, laser eye surgery or a general ocular examination.

	HIPP	KLINI	ARMO
ADULT			
	144 8001		



Nasal Specula

THUDICHUM

Thudichum Nasal Speculum are used to examine the nasal cavity. The blades are inserted into the nostrils to widen the nasal cavity to view the nasal structures.

	HIPP	KLINI	ARMO
NO.0			
18 x 8mm			A5400
NO.1			
20 x 10mm			A5401
NO.2			
22 x 12mm			A5402
NO.3			
24 x 13mm			A5403
NO.4			
26 x 13mm			A5404





19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au