

Vaginal Specula & Protoscopes

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Vaginal Specula

Grave	3
Kogan Endocervival Speculum	3
Pederson	3
Sims	3

Protoscopes

Kelly	4
-------	---

Uterine, Sounds

Sims Malleable	4
Sims Rigid	4



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

 <p>ANTON HIPPI Instruments & Implants</p>		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
 <p>KLINI™ instruments by anton hipp GERMANY</p>		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
 <p>ARMO surgical instruments</p>		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
 <p>Sayco first aid instruments</p>		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



Vaginal Specula



GRAVE

A Grave Vaginal Specula is a gynaecological instrument used for performing vaginal pelvic examinations.

	HIPP	KLINI	ARMO
SMALL			
20x75mm	5.010.01	K5.010.01	A5500
MEDIUM			
35x95mm	5.010.02	K5.010.02	A5501
LARGE			
35x115mm	5.010.03	K5.010.03	A5502
EXTRA LARGE			
38x128mm	5.010.04		




PEDERSON

Pederson Vaginal Specula is a gynaecological instrument used for performing vaginal pelvic examinations.

	HIPP	KLINI	ARMO
SMALL			
13x75mm	5.011.01	K5.011.01	A5503
MEDIUM			
22x100mm	5.011.02	K5.011.02	A5504
LARGE			
25x120mm	5.011.03	K5.011.03	A5505

Double Ended



Single Ended

SIMS

Sims Vaginal Specula is a gynaecological instrument. The Sims vaginal speculum is used to dilate a vagina so that the cervix and the vagina can be examined.

	HIPP	KLINI	ARMO
DOUBLE ENDED			
Small			A5511
Medium			A5512
Large			A5513
SINGLE ENDED			
Small			A5515
Medium			A5516
Large			A5517

KOGAN ENDOCERVIVAL SPECULUM

Kogan Endocervical Speculum is a gynaecological instrument used for viewing the distal endocervical canal. This single spring, locking instrument features opposing, fenestrated blades at the end of long jaws.

	HIPP	KLINI	ARMO
WITH SET SCREW & GAUGE			
3mm wide tips 24cm			A5506
WITH RATCHET, NO GAUGE			
3mm wide tips 28cm			A5507
7mm wide tips 28cm			A5508





Proctoscopes

KELLY

Kelly Proctoscope is used to examine the anal cavity, rectum or colon.

	HIPP	KLINI	ARMO
STANDARD			
21 x 50mm			A2390
MEDIUM			
27 x 50mm			A2391
LONG			
22 x 100mm			A2392




Uterine, Sounds

SIMS MALLEABLE

Sims Uterine Sound is an obstetric instrument used for probing and dilating the uterus as well as measuring the length of the cervical canal. The instrument has graduated centimetre markings.

	HIPP	KLINI	ARMO
Silver plated brass			A3819



SIMS RIGID

Sims Uterine Sound is an obstetric instrument used for probing and dilating the uterus as well as measuring the length of the cervical canal. The instrument has graduated centimetre markings.

	HIPP	KLINI	ARMO
Solid stainless steel			A3820

Instrument quality level



ANTON HIPF
Theatre (Germany)

KLINI
Germany

ARMO
Superior

Sayco
Ward



Vaginal Specula

05



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au