

Plaster Instruments

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Beeson Spreader	4
Bruns Plaster Scissors	3
Esmarch Plaster Scissors	3
Haglund-Stille Spreader	4
Hennig Spreader	4
Stille Cutting Shears	3
USA Model Spreader	4
Wolff Cast Breaker	3



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



WOLFF CAST BREAKER

Wolff Cast Breaker is used for breaking open plaster casts. It is a single spring, plier instrument. The jaws are curved and have horizontal serrations.

	HIPP	KLINI	ARMO
18cm	6.507.18		
24cm	6.507.24	K6.507.24	A6000



Smooth Jaw

Serrated Jaw

BRUNS PLASTER SCISSORS

Bruns Plaster Shears are heavy duty scissors/shears designed for use with plaster and fiberglass castings.

	HIPP	KLINI	ARMO
SMOOTH JAW			
24cm	1.559.23	K1.559.23	A6005
SERRATED JAW			
24cm	1.557.23	K1.557.23	A6006



ESMARCH PLASTER SCISSORS

Esmarch Plaster Shears are heavy duty angled scissors/shears designed for use with plaster and fiberglass castings.

	HIPP	KLINI	ARMO
20cm			A6007



STILLE CUTTING SHEARS

Stille Plaster Shears are large 2 handled shears used to remove heavy or large plaster casts. They have round handles to facilitate ease in the cast removal.

	HIPP	KLINI	ARMO
23cm			A6003
37cm			A6004



BEESON SPREADER

Beeson Plaster Spreader is used to aid in the removal of plaster casts. The blades of the spreader are closed by opening the handle and are inserted into a section of the cast. When squeezing the handle, the blades spread open thereby forcing the cast open.


	HIPP	KLINI	ARMO
30cm			A6002



HENNIG SPREADER

Hennig Plaster Spreader is used to aid in the removal of plaster casts. The blades of the spreader are closed by opening the handle and are inserted into a section of the cast. When squeezing the handle, the blades spread open thereby forcing the cast open.

	HIPP	KLINI	ARMO
27cm			A6001



USA MODEL SPREADER

USA Model Plaster Spreader is used to aid in the removal of plaster casts. The pointed blades of the spreader are closed by opening the handle and are inserted into a section of the cast. When squeezing the handle, the blades spread open thereby forcing the cast open.

	HIPP	KLINI	ARMO
23cm			A5998



HAGLUND-STILLE SPREADER

Haglund-Stille Plaster Spreader is used to aid in the removal of plaster casts. The tips of the spreader are closed by opening the handle and are inserted into a section of the cast. When squeezing the handle, the right angles tips spread open thereby forcing the cast open.

	HIPP	KLINI	ARMO
28cm			A5999

Instrument quality level



Theatre (Germany)

Germany

Superior

Ward



Plaster Instruments

05



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au