

Hooks Retractors & Elevators

BYDAND INSTRUMENT CATALOGUE
 VERSION 11



Hooks

Gillies	3
Hartman Cerumen	3
Joseph	3
Kilner	4
Spey (Snook Ovaryectomy)	3
Tracheal	4

Retractors

Delicate	6
Gelpi	5
Parsonnet Epicardial	5
Ragnell-Davis	6
Senn Miller	6
Weitlaner	5

Elevators

Septum	7
--------	---



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.



THEATRE (GERMANY)

Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.



GERMANY

Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.



SUPERIOR

The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.



WARD

Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.

Hooks

SPEY (SNOOK OVARECTOMY)

Spey hook is a long, thin, probe like instrument with a curved end and a buttoned tip, used in an ovariectomy procedure.

	HIPP	KLINI	ARMO
7cm shank			A3000
11cm shank			A3001

HARTMAN CERUMEN

Hartman Cerumen Hook allows for the easy extraction of hard wax and other obstacles from the ear canal.

	HIPP	KLINI	ARMO
16cm			A2999

GILLIES



Gillies Skin Hook is used to retract the skin during surgical procedures.

	HIPP	KLINI	ARMO
SHARP - 4MM HOOK			
18cm	6.152.02	K6.152.02	A3358
BLUNT - 4MM HOOK			
18cm	6.152.01		A3356

JOSEPH

Joseph Skin Hook is used to retract the skin during intranasal and pharyngeal surgical procedures. It has a round handle with dual sharp prongs on the end.

	HIPP	KLINI	ARMO
2MM			
16cm			A3351
5MM			
16cm			A3350





Hooks

KILNER

Kilner Skin Hook is primarily used during general surgeries when delicate tissue needs to be retracted.

	HIPP	KLINI	ARMO
SHARP - 2.5MM HOOK			
16cm	6.150.16	K6.150.16	A3357

TRACHEAL

Tracheal Hook is used in ENT surgery, usually to retract the trachea during a tracheotomy.

	HIPP	KLINI	ARMO
BLUNT - 1 PRONG			
16cm			A3364

Retractors



WEITLANER

Weitlaner Retractor is used for holding back tissue and exposing a surgical site. It has curved shanks that can be held open with a sprung ratchet mechanism.


	HIPP	KLINI	ARMO
SHARP - 2 X 3 TEETH			
10.5cm			A4159
SHARP - 3 X 4 TEETH			
14cm			A4157
BLUNT - 2 X 3 TEETH			
10.5cm			A4158
BLUNT - 3 X 4 TEETH			
14cm			A4156



GELPI

Gelpi Retractor is a ratcheted retractor used to hold back tissue during a surgical procedure. It has a single sharp prong on each tip that is turned outwards.

	HIPP	KLINI	ARMO
SHARP - 1 X 1 TEETH			
11cm	6.248.11	K6.248.11	A3353




PARSONNET EPICARDIAL

Parsonnet Epicardial Retractor is a small, sprung retractor used in coronary bypass surgery. It is available with 3X3 prongs and has thumb grips.

	HIPP	KLINI	ARMO
Mini	315 2506		
Small	315 3206		
Medium	315 3807		
Large	315 5011		





Sharp

2 prong

3 prong



Blunt


2 prong

3 prong



Sharp

Blunt



Retractors

DELICATE

Delicate Retractors are used to hold back tissue, skin and or muscle during surgery.

	HIPP	KLINI	ARMO
SHARP - 2 PRONG			
16cm			A3359
SHARP - 3 PRONG			
16cm			A3360
BLUNT - 2 PRONG			
16cm			A3361
BLUNT - 3 PRONG			
16cm			A3362

SENN MILLER

Senn Miller Retractor is a double-ended retractor used to retract surface tissue.

	HIPP	KLINI	ARMO
SHARP			
16cm	6.139.03	K6.139.03	A3352
BLUNT			
16cm	6.139.02	K6.139.02	A3355

RAGNELL-DAVIS

Ragnell-Davis Retractor is a double-ended retractor used to retract delicate surface tissue.

	HIPP	KLINI	ARMO
15cm			A3354

Retractors



ROUX

Roux Retractor is a double-ended handheld retractor with curved ends on each side. Used to hold back tissue during surgery.

	HIPP	KLINI	ARMO
14cm			A3376
15cm			A3377
16cm			A3378



BALFOUR ABDOMINAL

The Balfour Retractor is used to hold open the abdominal walls during surgery. The Balfour is made up of the frame with an arm on each end and a blade in the middle that can slide from side to side and be locked into place. The surgical retractor can either be self retaining or locked into place.

	HIPP	KLINI	ARMO
SHARP - 1 X 1 TEETH			
18cm spread		K6.277.18	



LANGENBECK

Langenbeck retractors are used to retract soft tissues and wound edges in order to visualize deeper structures during urological, plastic and general surgical procedures.

	HIPP	KLINI	ARMO
3x1.1cm			A3363

Elevators



Freer, Double ended

SEPTUM

Freer Septum Elevators are used in nasal septum surgeries.

	HIPP	KLINI	ARMO
BLUNT/SHARP			
18cm			A3370



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au