

Ring Cutters

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Standard
Walton

3
3



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

| | |
|--|---|
|  <p>ANTON HIPPI Instruments & Implants</p> | <p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p> |
|  <p>KLINI™ instruments by anton hipp GERMANY</p> | <p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p> |
|  <p>ARMO surgical instruments</p> | <p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p> |
|  <p>Sayco first aid instruments</p> | <p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p> |



Standard

RING CUTTER

Ring Cutters are plier gripped instruments used for cutting rings which have become difficult to remove. The blade is made from a hardened alloy and is completely enclosed and can quickly and safely remove rings from fingers. For use on soft metals only.

| | HIPP | KLINI | ARMO |
|------------------|------|-------|-------|
| 16cm | | | A5550 |
| 30mm Spare blade | | | A5560 |



Spare blade

Walton

WALTON RING CUTTER

Walton Ring Cutters are single gripped instruments used for cutting rings which have become difficult to remove. The blade is made from a hardened alloy and is completely enclosed and can quickly and safely remove rings from fingers. The safety lever slides protective a plate between finger and ring. For use on soft metals only.

| | HIPP | KLINI | ARMO |
|------------------|------|-------|-------|
| 16cm | | | A5555 |
| 20mm Spare blade | | | A5556 |



Spare blade



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au