

Dilators

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Trousseau	
Trousseau	3
Hegar	
Double Ended	3
Single Ended	3
Hawkin-Ambler	
Single Ended	4



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

 <p>ANTON HIPPI Instruments & Implants</p>	<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
 <p>KLINI™ instruments by anton hipp GERMANY</p>	<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
 <p>ARMO surgical instruments</p>	<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
 <p>Sayco first aid instruments</p>	<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



Trousseau

TROUSSEAU

Trousseau Dilator is a single spring instrument used to widen a tracheal incision to assist in tube insertion and placement.

	HIPP	KLINI	ARMO
14cm			A1023



Hegar

SINGLE ENDED

Hegar dilators are used to induce cervical dilation to gain entry to the uterus.

	HIPP	KLINI	ARMO
3mm			A1883
4mm			A1884
5mm			A1885
6mm			A1886
7mm			A1887
8mm			A1888
9mm			A1889
10mm			A1890
11mm			A1891
12mm			A1892
13mm			A1893
14mm			A1894
Set of 12			A1882



DOUBLE ENDED

Hegar dilators are used to induce cervical dilation to gain entry to the uterus.

	HIPP	KLINI	ARMO
1-2mm			A1901
3-4mm	5.090.04		A1903
5-6mm	5.090.06		A1905
7-8mm	5.090.08		A1907
9-10mm	5.090.10		A1909
11-12mm	5.090.12		A1911
13-14mm	5.090.14		A1913
15-16mm	5.090.16		A1915
17-18mm	5.090.18		A1917
19-20mm			A1919
21-22mm			A1921
23-24mm			A1923
25-26mm			A1925



Hawkin-Ambler

SINGLE ENDED

Hawkin-Ambler Dilators are single-ended uterine dilators with a rounded end and are used to facilitate entry into the uterus and to prevent uterine wall puncture.

	HIPP	KLINI	ARMO
3-6mm			A1005
4-7mm			A1006
5-8mm			A1007
6-9mm			A1008
7-10mm			A1009
8-11mm			A1010
9-12mm			A1011
10-13mm			A1012
11-14mm			A1013
12-15mm			A1014
13-16mm			A1015
14-17mm			A1016
15-18mm			A1017
16-19mm			A1018
17-20mm			A1019
18-21mm			A1020

Instrument quality level



Dilators

05



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au