

# Towel Clips

BYDAND INSTRUMENT CATALOGUE  
VERSION 11



Backhaus  
Cross Action  
Lorna

3  
3  
3



# FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p><b>THEATRE (GERMANY)</b></p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p><b>GERMANY</b></p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p><b>SUPERIOR</b></p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p><b>WARD</b></p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



### BACKHAUS

The Backhaus Towel Clip is a ratcheted instrument used to secure draping and towels which are used while a medical procedure is being conducted. It has curved tips which are designed to perforate the item they grip.

	HIPP	KLINI	ARMO
9cm			A5150
11cm			A5154
14cm			A5152



### CROSS ACTION

The Cross-Action Towel Clip is a sprung instrument used to secure draping and towels which are used while a medical procedure is being conducted. It has curved tips which are designed to perforate the item they grip.

	HIPP	KLINI	ARMO
8cm			5161



### LORNA

The Lorna Towel Clip is a ratcheted instrument used to secure draping and towels which are used while a medical procedure is being conducted. It has curved tips which are blunt.

	HIPP	KLINI	ARMO
14cm			A5170



19/9 Hudson Avenue  
Castle Hill NSW 2154  
Australia  
phone (02) 9659 9548  
email [admin@bydandmedical.com.au](mailto:admin@bydandmedical.com.au)

[bydandmedical.com.au](http://bydandmedical.com.au)