


Probes

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Applicators

Cotton (Standard)	6
Farrell Cotton	6
Jobson-Horne	6
Uckerman	6

Eye Probes

Eye Magnet & Foreign Body Loop	5
Eye Probes	4
Foreign Body Needle & Gouge	5
Foreign Body Needle & Spud	4
Foreign Body Spud & Gouge	5

Extractors

Schamberg	3
Unna with Teardrop ends	3
Unna with round ends	3

Lachrymal Probe

Lachrymal Probe	5
-----------------	---

Probes, Seekers & Directors

Alloy Handle	4
Stainless Steel	3

Splinter Probes

Disposable, Double Ended	4
--------------------------	---



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



Extractors

SCHAMBERG

Schamberg extractor is a dermal instrument used to extract blackheads.

	HIPP	KLINI	ARMO
10cm			A6027

UNNA WITH TEARDROP ENDS

Unna Extractor is a dermal instrument used to extract blackheads.

	HIPP	KLINI	ARMO
14cm			A6024

UNNA WITH ROUND ENDS

Unna Extractor is a dermal instrument used to extract blackheads.


	HIPP	KLINI	ARMO
14cm			A6025

Probes, Seekers & Directors

STAINLESS STEEL

Probes, seekers and Directors are long, slender instruments for exploring wounds or body cavities or passages.

	HIPP	KLINI	ARMO	SAYCO
BLUNT - 1mm Ø				
20cm	24.600.20			
BLUNT - 2mm Ø				
14cm			A2366	
20cm	24.602.20			
SHARP - 2mm Ø				
14cm			A2361	
DISSECTING - WITH EYE				
14cm				235
SHARP				
14cm				255
BLUNT				
14cm				253
BLUNT - ANGLED				
14cm				254



○
1mm
Diametre

○
2mm
Diametre



Probes, Seekers & Directors

ALLOY HANDLE

Probes, seekers and Directors are long, slender instruments for exploring wounds or body cavities or passages.

	HIPP	KLINI	ARMO	SAYCO
BLUNT				
14cm				237
BLUNT - ANGLED				
14cm				236
SHARP				
12cm				226



Splinter Probes

DISPOSABLE, DOUBLE ENDED

Splinter probes are used to push a splinter to the surface of the skin.

	HIPP	KLINI	ARMO	SAYCO
WITH CAPS				
11cm				519



Eye Probes

EYE PROBES

Eye probes are used to remove foreign bodies from the eye.

	HIPP	KLINI	ARMO	SAYCO
WITH CAPS				
11cm				520



FOREIGN BODY NEEDLE & SPUD

Foreign body needle and Spud are ophthalmic instruments commonly used to assist in the removal of foreign bodies in the eye.

	HIPP	KLINI	ARMO	SAYCO
			A6020	



Eye Probes

FOREIGN BODY NEEDLE & GOUGE

Foreign body needle and Spud are ophthalmic instruments commonly used to assist in the removal of foreign bodies in the eye.

	HIPP	KLINI	ARMO	SAYCO
WITH CAPS				
15cm			A6021	

FOREIGN BODY SPUD & GOUGE

Foreign body needle and Spud are ophthalmic instruments commonly used to assist in the removal of foreign bodies in the eye.

	HIPP	KLINI	ARMO	SAYCO
WITH CAPS				
13.5cm			A6022	

EYE MAGNET & FOREIGN BODY LOOP

Eye Magnet & Foreign Body Loop is used for removing objects from the eye. The magnet is designed to remove fragments of metal lodged near the surface of the eye. At the opposite end of the instrument is a foreign body loop used to encircle the object for extraction.


	HIPP	KLINI	ARMO	SAYCO
REVERSIBLE SCREW HANDLE				
11.5cm			A6023	

Lachrymal Probe

LACHRYMAL PROBE

The Lachrymal Probe Bowman is a very thin malleable rod with a rounded end used for clearing a nasolacrimal duct, or in procedures within the nasolacrimal duct, or in procedures within the lacrimal system.

	HIPP	KLINI	ARMO	SAYCO
00/0			A1780	
0/1			A1781	



Applicators

COTTON (STANDARD)

Cotton Standard Applicators have many uses in the medical environment. They can be used to clean the wound, as well skin surrounding it. They are also useful for gently swabbing a wound to obtain a sample for further testing that might be required to determine the presence of an infection. Cotton tipped applicators are useful in applying medicines or ointment to wounds as well.

	HIPP	KLINI	ARMO	SAYCO
20cm			A1024	



FARRELL COTTON

Farrell Cotton Applicator is commonly used in aural procedures.

	HIPP	KLINI	ARMO	SAYCO
14cm			A2380	



JOBSON-HORNE

Jobson-Horne Cotton Applicator is characterised by having a curette loop at one end and a flattened serrated tip at the other end. Commonly used in aural procedures.

	HIPP	KLINI	ARMO	SAYCO
12cm			A2381	
18cm			A2382	



UCKERMAN

Uckermann Cotton Applicator is characterised by having a curette loop at one end and a very fine spiral tip (<1mm) at the other end. Commonly used in aural procedures.

	HIPP	KLINI	ARMO	SAYCO
15cm			A2378	

Instrument quality level

 ANTON HIPPE Theatre (Germany)	 KLINI Germany	 ARMO Superior	 Sayco Ward
---	---	---	--



Probes

07



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au