

# Stethoscopes

BYDAND INSTRUMENT CATALOGUE  
 VERSION 11



Cardiology	2
Classic	2
Classic-S	2
Majestic	3
Pinnard Foetal	5
Spectrum, Lightweight Dual Head	5
Spectrum, Lightweight Single Head	4
Sprague Rappaport	4
Toonscope	3
Training	5



**CARDIOLOGY**

- Extra deep bell complete with non-chill ring on both bell and diaphragm sides.
- Sensitive diaphragm for low frequency sounds.
- Soft and flexible Y tubing with double leaf internal spring.
- Binaurals and chestpiece are made from stainless steel.
- Supplied with two pairs of soft deluxe eartips in small and large size.
- Spare diaphragm is included.

ABN	
<b>ADULT</b>	
Black	<b>SGSS-011-BK</b>
Burgundy	<b>SGSS-011-BD</b>
Hunter Green	<b>SGSS-011-HG</b>
Grey	<b>SGSS-011-GY</b>



**CLASSIC**

- Optimized acoustic response and rugged durability.
- Precisely machined chestpiece for outstanding performance.
- Double leaf internal binaurals for comfortable ear fit.
- Non-chill ring and diaphragm retaining ring for patient comfort.
- Extra set of soft sealing eartips and one spare diaphragm included.

ABN	
<b>ADULT</b>	
Black	<b>SGSS-021-BK</b>
Royal Blue	<b>SGSS-021-NB</b>
Burgundy	<b>SGSS-021-BD</b>
Hunter Green	<b>SGSS-021-HG</b>
Grey	<b>SGSS-021-GY</b>
Light Blue	<b>SGSS-021-LB</b>
<b>INFANT</b>	
Light Blue	<b>SGSS-023-LB</b>
<b>CHILD</b>	
Black	<b>SGSS-025-BK</b>



**CLASSIC-S**

- Excellent acoustic response and rugged durability.
- Precisely machined stainless steel chestpiece with outstanding performance.
- Tunable diaphragm for high and low frequencies.
- Non-chill ring for patient comfort
- Double internal leafspring binaurals for comfortable ear fit.
- Extra set of soft sealing eartips

ABN	
<b>ADULT</b>	
Black	<b>SGSS-021S-BK</b>
Royal Blue	<b>SGSS-021S-RB</b>
Burgundy	<b>SGSS-021S-BD</b>
Hunter Green	<b>SGSS-021S-HG</b>
Pink	<b>SGSS-021S-PK</b>



### MAJESTIC

- Fine acoustic response and a light weight design for comfort.
- Precision machined aluminum chestpiece.
- Double leaf internal spring binaurals for comfortable ear fit.
- Supplied with an extra set of soft sealing eartips and diaphragm.

ABN	
<b>ADULT</b>	
Black	<b>SGGS-071-BK</b>
Royal Blue	<b>SGGS-071-RB</b>
Burgundy	<b>SGGS-071-BD</b>
Hunter Green	<b>SGGS-071-HG</b>
Grey	<b>SGGS-071-GY</b>
Light Blue	<b>SGGS-071-LB</b>
Pink	<b>SGGS-071-PK</b>



### TOONSCOPE

- Single head stethoscope with seven fun cartoon faces in molded resin.
- Multi frequency design switches between low and high frequencies simply by altering pressure.
- Comes with extra pair of colour coordinated soft sealing eartips.
- Overall length is 770mm.


ABN	
<b>CHILD</b>	
Light Blue	<b>SGTS-181-LB</b>
Yellow	<b>SGTS-181-YL</b>
Pink	<b>SGTS-181-PK</b>
Light Green	<b>SGTS-181-LG</b>



### SPRAGUE RAPPAPORT

- Five interchangeable chestpiece fittings allow for custom user design.
- Heavy walled, double tube configuration for improved acoustic separation.
- Supplied with accessory puoch containing the bell chestpiece, 2 spare diaphragms, 1 pair mushroom eartips and 1 pair PVC eartips.

	ABN
<b>ADULT</b>	
Black	<b>SGRS-09-BK</b>
Navy Blue	<b>SGRS-09-NB</b>
Burgundy	<b>SGRS-09-BD</b>
Hunter Green	<b>SGRS-09-HG</b>
Grey	<b>SGRS-09-GY</b>
Pink	<b>SGRS-09-PK</b>
Purple	<b>SGRS-09-PU</b>
Red	<b>SGRS-09-RD</b>
Yellow	<b>SGRS-09-YL</b>
Light Blue	<b>SGRS-09-LB</b>



### SPECTRUM, LIGHTWEIGHT SINGLE HEAD

- Elegant style and lightweight comfort with affordable pricing.
- Colour selection allows for easy departmental coding or user preference.
- Lightweight chrome plated/polished aluminum chestpiece.
- Adjustable brass plated binuaurals with PVC eartips.

	ABN
<b>ADULT</b>	
Black	<b>SGLS-121-BK</b>
Royal Blue	<b>SGLS-121-RB</b>
Hunter Green	<b>SGLS-121-HG</b>
Grey	<b>SGLS-121-GY</b>
Pink	<b>SGLS-121-PK</b>
Purple	<b>SGLS-121-PU</b>
Red	<b>SGLS-121-RD</b>
Yellow	<b>SGLS-121-YL</b>



### SPECTRUM, LIGHTWEIGHT DUAL HEAD

- Elegant style and lightweight comfort with affordable pricing.
- Colour selection allows for easy departmental coding or user preference.
- Lightweight chrome plated/polished aluminum chestpiece.
- Adjustable brass plated binaurals with PVC eartips.

ABN	
<b>ADULT</b>	
Black	<b>SGLS-101-BK</b>
Royal Blue	<b>SGLS-101-RB</b>
Hunter Green	<b>SGLS-101-HG</b>
Grey	<b>SGLS-101-GY</b>
Pink	<b>SGLS-101-PK</b>
Purple	<b>SGLS-101-PU</b>
Red	<b>SGLS-101-RD</b>
Yellow	<b>SGLS-101-YL</b>
Magenta	<b>SGLS-101-MG</b>
Teal	<b>SGLS-101-TL</b>



### TRAINING

- Lightweight chrome plated/polished aluminum dual head chestpiece.
- Two sets of adjustable chrome plated binaurals.
- Two sets of black Y-tubing with Y-connector and extension tubing.

ABN	
<b>ADULT</b>	
Black	<b>SGLS-151-BK</b>



### PINNARD FOETAL

- Cone Shaped for superior acoustic response.
- 140mm High
- 50mm Diameter

ABN	
<b>ALUMINIUM</b>	
Silver	<b>SGPF-426</b>
<b>PLASTIC</b>	
Black	<b>SGPF-420</b>



19/9 Hudson Avenue  
Castle Hill NSW 2154  
Australia  
phone (02) 9659 9548  
email [admin@bydandmedical.com.au](mailto:admin@bydandmedical.com.au)

[bydandmedical.com.au](http://bydandmedical.com.au)