

Sphygmomanometers

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Aneroid

Econom - Economy Aneroid	2
Home Blood Pressure Kit	2

Digital

Rossmax BP Automatic Upper Arm Bluetooth Monitor	3
Rossmax BP Automatic Upper Arm Monitor (Fuzzy Logic)	3
Rossmax BP Automatic Wrist Monitor	2
Rossmax BP Large Cuff 34-46cm	3

Palm Held

Palm- HandHeld Aneroid	4
Palm Lite - Hand Held Aneroid	4

Regal

Desk - Aneroid Sphygmomanometer Regal	4
Mobile - Aneroid Sphygmomanometer Regal	5
Wall Mounted - Aneroid Sphygmomanometer Regal	5



Aneroid

ECONOM - ECONOMY ANEROID

- Economically priced precision sphygmomanometer
- 300mmHg manometer with no pin stop mechanism.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
 - Supplied with zippered leatherette storage case.

	ABN
Kit	SGAS-051



HOME BLOOD PRESSURE KIT

Precision sphygmomanometer for professional and home use

- Designed for professional clinical use or for easy home self monitoring by individuals.
- Comes with Nylon D-ring cuff and attached stethoscope.
- D-ring cuff fits arms 25 to 35cm in circumference.
- 300mmHg manometer with no pin stop mechanism.
- Complete with heavy duty 2 tube inflation bladder, bulb and valve.
- Supplied with zippered leatherette storage case.

	ABN
Kit	SGAS-061

Digital

ROSSMAX BP AUTOMATIC WRIST MONITOR

Features:

- Smart Sense technology
- Irregular Heartbeat Detection (IHB)
- Hypertension risk indication
- Average of last 3 readings
- 120 Memories for 2-users with date and time indication
- Movement detection
- Jumbo size screen
- PC link
- Lifetime calibration
- One-touch automatic operation



	ROSSMAX
Kit	RMBQ75

Digital



ROSSMAX BP AUTOMATIC UPPER ARM MONITOR (FUZZY LOGIC)

Automatic BP monitor with lifetime calibration.

- Real Fuzzy Logic Technology
- Cuff Detection
- Movement Detection
- Hypertension Risk Indication
- Irregular Heartbeat Detection
- Average of Last 3 readings
- 120 Memories for 2 users with date and time
- Large LCD
- One Touch Operation

ROSSMAX

Kit **RMX3**



ROSSMAX BP AUTOMATIC UPPER ARM BLUETOOTH MONITOR

Automatic BP Monitor with Arrhythmia detection and lifetime calibration.

- Bluetooth low energy technology
- Arrhythmia Detection (ARR)
- Atrial Fibrillation Detection
- Premature Contraction Detection
- Movement Detection
- Cuff Detection
- Hypertension Risk Indication
- Average of last 3 readings
- 7 day AM/PM Averaging
- Patented Latex Free Universal Cone Cuff 24-40cm

ROSSMAX

Kit **RMX5BT**



ROSSMAX BP LARGE CUFF 34-46CM

No more annoying rotating and sliding when wrapping the cuff! The ergonomic, cone-shaped design of Rossmax Cone cuffs promotes blood pressure measurement to a whole new easy and delightful experience! This cuff is Latex free and Easy to apply and is available in small, medium and large.

ROSSMAX

Kit **BPCF001**



Palm Held

PALM- HANDHELD ANEROID

- Deluxe handheld model sphygmomanometer
- Deluxe Sphygmomanometer palm style with 300mmHg with no pin stop mechanism.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 1 tube inflation bladder.
 - Chrome plated inflation valve.
 - Unit supplied with deluxe zippered nylon storage case.

	ABN
Kit	SGAS-071LT



PALM LITE - HAND HELD ANEROID

- Handheld sphygmomanometer with precision performance and a modern lightweight design
- Precision crafted 300mmHg palm style manometer with enhanced shock resistance housing.
 - Quick release air valve button.
 - Quick release luer slip connector.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 1 tube inflation bladder.
 - Unit supplied with deluxe zippered nylon storage case.

	ABN
Kit	SGAS-071



Regal

DESK - ANEROID SPHYGMOMANOMETER REGAL

- Portable clock type Sphygmomanometer for professionals.
- Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
 - Easy to read black background with white writing.
 - Special non-reflective manometer glass.
 - Extra strong reinforced base design.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
 - Convenient under manometer cuff storage.
 - 120cm extendable coiled tubing included.

	REGAL
Kit	SGAS-081D

Regal

MOBILE - ANEROID SPHYGMOMANOMETER REGAL

- Clock type Sphygmomanometer on mobile Floor Stand.
- Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
 - Easy to read black background with white writing.
 - Special non-reflective manometer glass.
 - Two-piece height adjustable stand.
 - Five caster wheel base ensures easy movement and high stability.
 - Built-in cuff storage.
 - Stand is finished in durable chrome and clinical white enamel.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
 - 240cm extendable coiled tubing included.
 - All necessary assembly tools included.

REGAL

Kit **SGAS-081M**



WALL MOUNTED - ANEROID SPHYGMOMANOMETER REGAL

- Portable clock type Sphygmomanometer for professionals.
- Alternative presentation of the CLOCK aneroid sphygmomanometer.
 - Attaches to 180° swivel wall mounting bracket.
 - Large 150mm faceplate 300mmHg manometer with no pin stop mechanism.
 - Easy to read black background with white writing.
 - Special non-reflective manometer glass.
 - Extra strong reinforced storage basket.
 - Includes Adult size nylon fabric cuff with velcro.
 - Complete with heavy duty 2 tube inflation bladder, bulb and valve.
 - 240cm extendable coiled tubing included.

REGAL

Kit **SGAS-081W**





19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au