

Laryngoscopes

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Standard Laryngoscope Sets

Handles	3
McIntosh Blades	4
McIntosh Set (Standard)	3
Miller Blades	4
Miller Set (Standard)	3
Spare Standard Bulb	4

Fibre Optic Laryngoscope Sets

Fibre Optic McIntosh Blades	6
Fibre Optic McIntosh Set	5
Fibre Optic Miller Blades	6
Fibre Optic Miller Set	5
Handles	5
Spare Fibre Optic Light Bulb	6



FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

		<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
		<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
		<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
		<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



Standard Laryngoscope Sets



MCINTOSH SET (STANDARD)

Laryngoscope features

- Stainless steel (416L)
- Laryngoscope Sets are supplied with blades and handle in protective hard case
- McIntosh Set includes McIntosh blades 1 to 4 and a 32mm handle
- Sets are supplied excluding batteries

	ARMO
Set	A0204



MILLER SET (STANDARD)

Laryngoscope features

- Stainless steel (416L)
- Laryngoscope Sets are supplied with blades and handle in protective hard case
- Miller Set includes Miller blades 0 to 3 and a 32mm handle
- Sets are supplied excluding batteries

	ARMO
Set	A0614

19mm Ø
Handle



32mm Ø
Handle





HANDLES

Handle features

- Handles are interchangeable with McIntosh and Miller blades (all sizes)
- 19mm and 32mm diameter handles with sure grip surfacing
- Compatible with most standard hook-on blades that conform to ASTM(F965)

	ARMO
19mm Ø Handle	A0195
32mm Ø Handle	A0192



Standard Laryngoscope Sets

MCINTOSH BLADES

- Blade features
- Surgical grade stainless steel
 - Autoclavable to 134° C after removal of bulb
 - Compatible with most standard hook-on handles that conform to ASTM(F965)

	ARMO
Blade No #0	A0191
Blade No #1	A0187
Blade No #2	A0188
Blade No #3	A0189
Blade No #4	A0190

MILLER BLADES

- Blade features
- Surgical grade stainless steel
 - Autoclavable to 134° C after removal of bulb
 - Compatible with most standard hook-on handles that conform to ASTM(F965)

	ARMO
Blade No #0	A0615
Blade No #1	A0616
Blade No #2	A0617
Blade No #3	A0618

SPARE STANDARD BULB

	ARMO
Spare bulb for F.O. Blades #0	A0223
Spare bulb for F.O. Blades #1-5	A0194



Fibre Optic Laryngoscope Sets



FIBRE OPTIC MCINTOSH SET

Laryngoscope features

- Laryngoscope Sets are supplied with blades and handle in protective hard case.
- McIntosh Set includes McIntosh blades 1 to 4 and a 32mm handle
- Sets and handles are supplied without batteries

	ARMO
Set	A0201




FIBRE OPTIC MILLER SET

Laryngoscope features

- Stainless steel (416L)
- Laryngoscope Sets are supplied with blades and handle in protective hard case
- Miller Set includes Miller blades 0 to 3 and a 32mm handle
- Sets are supplied excluding batteries

	ARMO
Set	A0202




HANDLES

Handle features

- Handles are interchangeable with McIntosh and Miller blades (all sizes)
- 32mm diameter handle with sure grip surfacing
- Handles are compatible with most F.O. hook-on blades conforming to the ISO 7376 (green standard)

	ARMO
F.O. Handle	A0203



Fibre Optic Laryngoscope Sets

FIBRE OPTIC MCINTOSH BLADES

Blade features

- Surgical grade stainless steel
- Autoclavable to 134° C
- Handles are compatible with most F.O. hook-on blades conforming to the ISO 7376 (green standard)

	ARMO
F.O. McIntosh Blade #1	A0716
F.O. McIntosh Blade #2	A0717
F.O. McIntosh Blade #3	A0718
F.O. McIntosh Blade #4	A0719



FIBRE OPTIC MILLER BLADES

Blade features

- Surgical grade stainless steel
- Autoclavable to 134° C
- Handles are compatible with most F.O. hook-on blades conforming to the ISO 7376 (green standard)

	ARMO
F.O. Miller Blade #0	A0815
F.O. Miller Blade #1	A0816
F.O. Miller Blade #2	A0817
F.O. Miller Blade #3	A0818



SPARE FIBRE OPTIC LIGHT BULB

	ARMO
	A0226

Instrument quality level



Theatre (Germany)



Germany



Superior



Ward



Laryngoscopes

07



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au