

## Tuning Forks, Reflex Hammers, Penlights, Tapes and Miscellaneous

BYDAND INSTRUMENT CATALOGUE  
VERSION 11



<b>Diagnostic Penlights</b>	
Deluxe Penlight	3
Economy Penlight	3
<b>Measuring Tapes</b>	3
<b>Pin Wheel</b>	
Wartenburg Pin Wheel	4
<b>Reflex Hammers</b>	
Buck (with needle & brush)	4
National/London Hospital	4
Taylor	4
<b>Tuning Forks</b>	
Armo Alloy	5
Armo Alloy Boxed SET	5
Hartmann Stainless	5
Lucae	6
Rydel-Seiffer	5



# FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

<b>ANTON HIPPI</b> Instruments & Implants	<p><b>THEATRE (GERMANY)</b></p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
<b>KLINI</b> instruments by anton hipp GERMANY	<p><b>GERMANY</b></p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
<b>ARMO</b> surgical instruments	<p><b>SUPERIOR</b></p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
<b>Sayco</b> first aid instruments	<p><b>WARD</b></p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>



## Diagnostic Penlights

### DELUXE PENLIGHT

- High intensity light beam
- Modern look chrome finish
- Pocket clip on/off press switch
- Pupil gauge and slide rule
- Includes 2 x AAA batteries

	<b>ABN</b>
Chrome finish	<b>SGPL-900-212</b>

### ECONOMY PENLIGHT

- High intensity light beam
- Stylish white finish with pupil gauge
- Pocket clip on/off press switch
- Includes 2 x AAA batteries

	<b>ABN</b>
Plastic	<b>SGPL-900-213</b>

## Measuring Tapes

### MEASURING TAPE

- Stylish white compact housing
- Push button tape rewind
- Measurement in inches and centimetres (60"/150cm)

	<b>ABN</b>
Plastic	<b>SGTM-900-300</b>





Instrument quality level

ANTON HIPPI  
Theatre (Germany)

KLINI  
Germany

ARMO  
Superior

SAYCO  
Ward



## Pin Wheel

### WARTENBURG PIN WHEEL

The Wartenburg Pin Wheel is a neurological testing instrument used for testing nerve response when pressure is applied to the skin.

	HIPP	KLINI	ARMO
18cm			A2009



## Reflex Hammers

### NATIONAL/LONDON HOSPITAL

	HIPP	KLINI	ARMO	SAYCO
38cm				232



### BUCK (WITH NEEDLE & BRUSH)

	HIPP	KLINI	ARMO	SAYCO
18cm				244



### TAYLOR

	HIPP	KLINI	ARMO	SAYCO
20cm				243

Instrument quality level



## Tuning Forks

### ARMO ALLOY

	HIPP	KLINI	ARMO
<b>STANDARD WITH WEIGHTS</b>			
Standard 128Hz			A1538
Standard 256Hz			A1539
<b>STANDARD</b>			
Standard 128Hz			A1538
Standard 256Hz			A1539
Standard 2048Hz			A1542

### ARMO ALLOY BOXED SET

	HIPP	KLINI	ARMO
Set of 5			A1537

### HARTMANN STAINLESS

	HIPP	KLINI	ARMO
<b>STANDARD</b>			
Standard 128Hz			A1530
Standard 128Hz (A)			A1531
<b>STANDARD</b>			
Standard 512Hz			A1532
Standard 1024Hz			A1533

### RYDEL-SEIFFER

	HIPP	KLINI	ARMO
64 - 128Hz			A1543





Instrument quality level



## Tuning Forks

### LUCAE

	HIPP	KLINI	ARMO
<b>STANDARD WITH FOOT</b>			
Standard 128Hz			A1550
Standard 256Hz			A1551
Standard 512Hz			A1552
Standard 1024Hz			A1553
Standard 2048Hz			A1554

Instrument quality level



Theatre (Germany)

Germany

Superior

Ward



## Tuning Forks, Reflex Hammers, Penlights, Tapes and Miscellaneous

07



19/9 Hudson Avenue  
Castle Hill NSW 2154  
Australia  
phone (02) 9659 9548  
email [admin@bydandmedical.com.au](mailto:admin@bydandmedical.com.au)

[bydandmedical.com.au](http://bydandmedical.com.au)