




# Single Use Instruments

BYDAND INSTRUMENT CATALOGUE  
VERSION 11



Contents (by kind)	2
Artery Forceps	7
Dressing Forceps	5
Instrument Packs	10
Needle Holders	3
Scalpel Handle	9
Scissors	4
Splinter Forceps	5
Sponge Holding Forceps	7
Suture Packs	9
Tissue Forceps	6
Uterine Forceps	8
Uterine Sound	8




## SINGLE USE INSTRUMENTS

Bydand Medical offer a comprehensive range of quality instruments in a single use presentation for those medical professionals that prefer the convenience of using instruments that are already sterilised.

Our **Sayco Sterile** range is a much welcomed addition to compliment our existing instrument ranges.

### SAYCO STERILE INSTRUMENTS ARE:

- Sterile, ready-to-use instruments.
- Supplied in easy-to-open durable clear plastic pouches.
- Blue colour coded handles for instant identification as single use.
- Available in packs of 10.
- **GUARANTEED QUALITY.** As with all the instruments that bear our brand name we guarantee that our product will be free of functional defects in workmanship and materials when used correctly in its intended purpose.



## Contents

### Artery Forceps

Crile	8
Halstead Mosquito - Micro	7
Halstead Mosquito	7
Hartman Alligator	8
Kelly	7

### Dressing Forceps

Adson	5
Dressing Forceps	5
McIndoe	5

### Instrument Packs

Dressing	10
Gillies	10
Semken	10

### Needle Holders

Crile-Wood	3
Derf	3
Mayo-Hegar	3

### Scalpel Handle

Scalpel Handle #3	9
-------------------	---

### Scissors

Dressing	4
Iris	4
Sims Uterine	4

### Splinter Forceps

Hunter	5
--------	---

### Sponge Holding Forceps

Ramplcy	7
---------	---

### Suture Packs

Micro Suture Pack	9
Suture Pack with SH/BL Scissors	9
Suture Pack with SH/SH Scissors	10

### Tissue Forceps

Adson	6
Adson Micro	6
Gillies	6
Semken	6

### Uterine Forceps

Teale Vulsellum	8
-----------------	---

### Uterine Sound

Sims	8
------	---



## Needle Holders

### MAYO-HEGAR

Mayo-Hegar Needle Holders are amongst the most popular locking needle holders used. They have a cross-hatched pattern on the jaws which provides a secure grip on suture needles. Our Single Use range are 15cm in length and are sold in packs of 10.

	SAYCO	Pack
12.5cm	BSU401	10

### DERF

Derf Needle Holders feature short jaws and a fine cross hatched pattern on the jaws that are suitable for small suture needles (commonly 5-0 and 6-0 sutures). Often utilized in ophthalmic, dental, or plastic surgery procedures. Our Single Use range are 12.5cm in length and are sold in packs of 10.

	SAYCO	Pack
12.5cm	BSU400	10

### CRILE-WOOD

Crile-Wood Needle Holders are ideally suited to small to medium size needles. Our Single Use range are 15cm in length and are sold in packs of 10.

	SAYCO	Pack
15cm	BSU405	10





## Scissors



### IRIS

Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications. Our Single Use range are available in straight and curved in an 11.5cm length - they are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
11.5cm	<b>BSU300</b>	10
<b>CURVED</b>		
11.5cm	<b>BSU301</b>	10



### DRESSING

Surgical or operating scissors are used for general cutting and dissecting of tissue, bandages and gauze. Our Single Use range is available in 10cm length and packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
12.5cm sharp/blunt	<b>BSU304</b>	10



### SIMS UTERINE

Sims (Uterine) scissors have long blades, primarily used for deep tissue dissection. These scissors are available with blunt tips with a straight blade in a 23cm length. They are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
23cm	<b>BSU321</b>	10



## Dressing Forceps

### DRESSING FORCEPS

Dressing Forceps are designed to hold gauze and other dressings. They are for general use including grasping dressings, suture material, tissue or foreign bodies in wounds. Our Single Use range are 12cm in length and are sold in packs of 10.

	SAYCO	Pack
12cm	BSU101	10

### ADSON

Adson Dressing Forceps are used for holding materials such as cotton and gauze, changing dressings, or packing wounds. They have a wide thumb grasp for improved control. Our Single Use range are 12cm in length and are sold in packs of 10.

	SAYCO	Pack
12cm	BSU104	10

### MCINDOE

McIndoe Forceps are used for holding materials such as cotton and gauze, changing dressings, or packing wounds. These forceps have a stop peg to ensure that when used for holding tissue the grip is not crushing. Our Single Use range are 15cm in length and are sold in packs of 10.

	SAYCO	Pack
15cm	BSU107	10





## Splinter Forceps

### HUNTER

Hunter Forceps are characterised by having an arrowhead pointed serrated tip and are used for the removal of splinters and other foreign bodies or to hold and manoeuvre dressings or soft tissue during procedures. Our Single Use range are 11.5cm in length and are sold in packs of 10.

	SAYCO	Pack
11.5cm	BSU108	10





## Tissue Forceps

### ADSON

Adson Tissue Forceps are used for holding and manipulating delicate tissue. Our Single Use range are 12cm in length and the jaw has 1x2 teeth - they are sold in packs of 10.

	SAYCO	Pack
<b>1X2 TEETH</b>		
12cm	<b>BSU103</b>	10

### ADSON MICRO

Adson Tissue Forceps are used for holding and manipulating delicate tissue. Our Single Use range are 12cm in length and the jaw has 1x2 teeth - they are sold in packs of 10.

	SAYCO	Pack
<b>1X2 TEETH</b>		
12cm	<b>BSU105</b>	10

### GILLIES

Gillies Tissue Forceps are used in small surgical sites to grasp tissue. They have fine, narrow jaws that are commonly used in plastic surgery. Our Single Use range are 15cm in length and the jaw has 1x2 teeth - they are sold in packs of 10.

	SAYCO	Pack
<b>1X2 TEETH</b>		
15cm	<b>BSU106</b>	10

### SEMKEN

Semken Tissue Forceps are lightweight, narrow forceps used for gripping fine membranes or single layer tissues.

	SAYCO	Pack
<b>1X2 TEETH</b>		
13cm	<b>BSU109</b>	10



## Sponge Holding Forceps

### RAMPLEY

Rampley Sponge Forceps are used to hold sponges and swabs. The jaws are elongated and contain serrated tips to hold swabs or sponges. Our Single Use range are 25cm in length - they are sold in packs of 10.

	SAYCO	Pack
25cm	BSU501	10



## Artery Forceps

### HALSTEAD MOSQUITO - MICRO

Halsted Micro Mosquito Forceps are ratcheted forceps and are used as a haemostat for small blood vessels. Our Single Use range are 12.5cm in length and available in straight and curved- they are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
12.5cm	BSU201	10
<b>CURVED</b>		
12.5cm	BSU202	10



### HALSTEAD MOSQUITO

Halsted Mosquito Forceps are ratcheted forceps and are used as a haemostat for small blood vessels. Our Single Use range are 12.5cm in length and available in straight and curved- they are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
12.5cm	BSU203	10
<b>CURVED</b>		
12.5cm	BSU204	10



### KELLY

Kelly Artery Forceps are ratcheted forceps used for clamping large blood vessels or manipulating heavy tissue. Our Single Use range are 14cm in length and available in straight and curved- they are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
14cm	BSU206	10
<b>CURVED</b>		
14cm	BSU205	10





## Artery Forceps

### CRILE

Crile Artery Forceps are ratcheted forceps used for clamping blood vessels or tissue before cauterization or ligation. Our Single Use range are 14cm in length straight - they are sold in packs of 10.


	SAYCO	Pack
<b>STRAIGHT</b>		
14cm	BSU207	10



### HARTMAN ALLIGATOR

The Hartman Alligator Forceps are medical forceps for removing foreign bodies and our Single Use range are 8.5cm in length and are sold in packs of 10.

	SAYCO	Pack
<b>STRAIGHT</b>		
8.5cm	BSU208	10



## Uterine Forceps

### TEALE VULSELLUM

Teale Vulsellum Forceps are curved ratcheted forceps used in obstetric and gynaecological surgery. Our Single Use range are 25cm in length - they are sold in packs of 10.

	SAYCO	Pack
<b>CURVED</b>		
25cm	BSU102	10



## Uterine Sound

### SIMS

Sims Uterine Sound is an obstetric instrument used for probing and dilating the uterus as well as measuring the length of the cervical canal. The instrument has graduated centimetre markings. Our Single Use range are 28cm in length - they are sold in packs of 10.

	SAYCO	Pack
<b>RIGID</b>		
28cm	BSU061	10





## Scalpel Handle

### SCALPEL HANDLE #3

Our Single Use Scalpel Handle # 3 are sold in packs of 10

	SAYCO	Pack
#3	BSU009	10



## Suture Packs

### MICRO SUTURE PACK

Included in this pack is

- 1x Iris Scissors Straight 11.5cm
- 1x Adson Tissue Forceps 1x2 Teeth 12cm
- 1x Derf Needle Holder 12.5cm

They are sold in packs of 10.

	SAYCO	Pack
Micro Suture Pack	BSU901	10



### SUTURE PACK WITH SH/BL SCISSORS

Included in this pack is

- 1x Wagner Scissors SH/BL Straight 12.5cm
- 1x Adson Tissue Forcep 1x2 Teeth 12cm
- 1x Mayo Hegar Needle Holder 15cm

They are sold in packs of 10.

	SAYCO	Pack
Suture Pack (SH/BL Scissors)	BSU902	10





## Suture Packs

### SUTURE PACK WITH SH/SH SCISSORS

Included in this pack

- 1x Wagner Scissors SH/SH Straight 12.5cm
- 1x Adson Tissue Forcep 1x2 Teeth 12cm
- 1x Mayo Hegar Needle Holder 15cm

They are sold in packs of 10.

	SAYCO	Pack
Suture Pack (SH/SH Scissors)	BSU903	10



## Instrument Packs

### GILLIES

Included in this pack

- 1x Iris Scissors - straight, 11.5cm
- 1x Gillies Tissue forcep 1x2
- 1x Mayo Hegar Needle Holder, 15cm

They are sold in packs of 10.

	SAYCO	Pack
Gillies Pack	BSU904	10



### SEMKEN

Included in this pack

- 1x Iris Scissor - straight, 11.5cm
- 1x Semken Tissue Forceps 1x2, 12cm

They are sold in packs of 10.

	SAYCO	Pack
Semken Pack	BSU905	10



## Instrument Packs



### DRESSING

Included in this pack

- 1 x Scissors SH/BL - straight, 12.5cm
  - 1 x Standard Dressing - blunt, 12cm
  - 1 x Mayo Hegar Needle Holder, 15cm
- They are sold in packs of 10.

	SAYCO	Pack
Dressing Pack	BSU906	10





19/9 Hudson Avenue  
Castle Hill NSW 2154  
Australia  
phone (02) 9659 9548  
email [admin@bydandmedical.com.au](mailto:admin@bydandmedical.com.au)

[bydandmedical.com.au](http://bydandmedical.com.au)