
SURGICRYL 910 POLYGLACTINE



Туре	Braided and coated multifilament				
Composition	Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide				
Coating	Poly-glycolide-co-L-lactide and calcium stearate (<1%)				
Colour	Violet				
Tissue reaction	Minimal				
Absorption	The hydrolytic action by which the material is broken down results in total absorption between 56 and 70 days. Between 40% and 50% of tensile strength remain after 21 days.				
Presentation	Box with 12 sutures				

USP 310 EP 2 REF 15201519

SURGICA VI. 910

Strate Similaria Aboutbable Subme 2/0 mirela nectorma.

For quality & safety





CHARACTERISTICS

High tensile strength

 $\overline{\mathbf{V}}$ Synthetic

Breaks down by hydrolysis

 $\sqrt{}$ Braided structure

Coated multifilament

Hermetically sealed packing

BENEFITS

Excellent tissue support

Minimal tissue reaction

Predictable absorption rate

Excellent handling properties

High flexibility High tensile strength Excellent knot security

Soft passage through the tissues with minimal sawing,

tissue drag and trauma

Good knot tie down / adjustability

Reduced capillary action

Guaranteed seal and product sterility

General soft tissue closing **INDICATIONS**

and/or ligation

especially in General surgery

Skin closure

Gastrointestinal surgery

Plastic surgery

Ophthalmic surgery

Gynaecology Obstetrics Orthopaedics

Urology

CONTRAINDICATIONS Cardiovascular surgery


Neurosurgery

Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised

immune system).

TENSILE STRENGTH



ABSORPTION



OPHTHALMIC SURGERY



Side cutting double arm			Metric	USP	Length		Order number
1/2 Circle	2 x HSP 8 mm 200 µm		0.4	8/0	45 cm	violet	O910081
	2 x HSP 8 mm 200 μm 300 μm		0.5	7/0	45 cm	violet	O91007
			0.7	6/0	45 cm	violet	O910061
	2 x HSP 8 mm 300 μm		0.7	6/0	45 cm	undyed	O910061U
Side cutting	Metric	USP	Length		Order number		
3/8 Circle	2 x DSP 6 mm 200 μm 220 μm		0.4	8/0	45 cm	violet	O91008
			0.5	7/0	45 cm	violet	O910072
	2 x DSP 11 mm 350 μm	\supset	1	5/0	45 cm	violet	O910051
Side cutting spatula double armed Metric USP Length Order num						Order number	
1/4 circle	2 x VSP 8 mm 300 μm		0.7	6/0	45 cm	violet	O91006
		0.7	6/0	20 cm	undyed	O91006U-20	
	•		0.7	6/0	45 cm	undyed	O91006U
	2 x VSP 8 mm		1	5/0	45 cm	violet 	O91005
	350 µm		1	5/0	45 cm	undyed	O91005U

