

Vet Surgical Kits

BYDAND INSTRUMENT CATALOGUE
VERSION 11



Canine Spey Kits	3
Feline Spey Kits	4
Ophthalmic kits	7
Soft Tissue Kits	6
Suturing Kits	5




FOUR CHOICES OF QUALITY

Our selection of HIPPI, KLINI, ARMO, and SAYCO instruments offer unparalleled flexibility, quality and value. Our multi-tiered instrument quality selection recognises and meets a variety of clinical needs and budgets.

  VINZENZ SATTLER	<p>THEATRE (GERMANY)</p> <p>Our state of the art range of "Theatre Quality" Instruments. All Hipp instruments are manufactured to exacting standards at the fully accredited Anton Hipp GMBH facility located in southern Germany. Crafted from the finest German stainless steel forgings for the utmost reliability, durability, consistency and precision. Incorporating the very latest technological innovations, most of these instruments are available with Tungsten Carbide inserts for added precision and lifespan.</p>
	<p>GERMANY</p> <p>Our new range of German theatre quality medical instruments. KLINI instruments are manufactured for us under licence by Anton Hipp. The HIPPI name is synonymous with quality that has been proven time and time again with over 40 years experience in the medical instrument industry. Our new KLINI range offers German quality at a competitive price.</p>
	<p>SUPERIOR</p> <p>The workhorse of our three quality ranges. This accredited range was borne out of the need for high quality dependable instrumentation in today's cost-sensitive environment. ARMO Instruments are crafted from the finest quality forgings and made to our specification in state of the art fully accredited manufacturing facilities in Pakistan. Most ARMO instruments are available in a smooth satin low reflective finish.</p>
	<p>WARD</p> <p>Our standard "Ward Quality" line. Economically priced and highly popular. This range is intended for use in applications where surgical precision is not the primary selection criterion.</p>




Canine Spey Kits



CANINE SPEY KITS			A9010-G1		A9010-G2		A9010-G3	
<i>Instrument Inclusions</i>			<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>
Rochester Carmalt Artery Forcep Straight	16cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Crile Artery Forceps Straight	14cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Crile Artery Forceps Curved	14cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Halstead-Mosquito Artery Forceps Curved	12.5cm		ARMO	4 x	ARMO	4 x	ARMO	4 x
Adson Dressing Forcep Standard	12cm		KLINI	1 x	KLINI	1 x	ARMO	1 x
Adson Brown Tissue Forcep	12cm		KLINI	1 x	KLINI	1 x	ARMO	1 x
Allis Tissue Forceps 4x5 teeth	15cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Spey Hook (Snook Ovaryectomy) shank	11cm		ARMO	1 x	ARMO	1 x	ARMO	1 x
Mayo Scissors Straight	14.5cm		HIPP	1 x	KLINI	1 x	ARMO	1 x
Metzenbaum Scissors Curved blunt/blunt	15.5cm		HIPP	1 x	KLINI	1 x	ARMO	1 x
Scalpel Handle (standard grip)	no.4		ARMO	1 x	ARMO	1 x	ARMO	1 x
Backhaus Towel Clips	9cm		ARMO	4 x	ARMO	4 x	ARMO	4 x



Feline Spey Kits



FELINE SPEY KITS			A9013-G1		A9013-G2		A9013-G3	
<i>Instrument Inclusions</i>			<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>
Rochester Carmalt Artery Forcep Straight	16cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Crile Artery Forceps Straight	14cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Crile Artery Forceps Curved	14cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Halstead-Mosquito Artery Forceps Curved	12.5cm		ARMO	4 x	ARMO	4 x	ARMO	4 x
Adson Dressing Forcep Standard	12cm		KLINI	1 x	KLINI	1 x	ARMO	1 x
Adson Brown Tissue Forcep	12cm		KLINI	1 x	KLINI	1 x	ARMO	1 x
Allis Tissue Forceps 4x5 teeth	15cm		ARMO	2 x	ARMO	2 x	ARMO	2 x
Spey Hook (Snook Ovarectomy) shank	7cm		ARMO	1 x	ARMO	1 x	ARMO	1 x
Mayo Scissors Straight	14.5cm		HIPP	1 x	KLINI	1 x	ARMO	1 x
Metzenbaum Scissors Curved blunt/blunt	15.5cm		HIPP	1 x	KLINI	1 x	ARMO	1 x
Scalpel Handle standard grip	no.4		ARMO	1 x	ARMO	1 x	ARMO	1 x
Backhaus Towel Clips	9cm		ARMO	4 x	ARMO	4 x	ARMO	4 x




Suturing Kits



CANINE SPEY KITS			A9016-G2		A9016-G3	
Instrument Inclusions			Brand	Qty	Brand	Qty
Crile Artery Forceps Curved	14cm		ARMO	2 x	ARMO	2 x
Adson Dressing Forcep Standard	12cm		KLINI	1 x	ARMO	1 x
Tissue Forcep 1x2 Teeth	14cm		ARMO	1 x	ARMO	1 x
Mayo Scissors Straight	14.5cm		KLINI	1 x	ARMO	1 x
Scalpel Handle standard grip	no. 4		ARMO	1 x	ARMO	1 x
Backhaus Towel Clips	9cm		ARMO	4 x	ARMO	4 x




Soft Tissue Kits



SOFT TISSUE KITS			A9018-G1		A9018-G2		A9018-G3	
<i>Instrument Inclusions</i>			<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>	<i>Brand</i>	<i>Qty</i>
Rochester Carmalt Artery Forcep Straight	16cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Rochester Carmalt Artery Forcep Curved	16cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Crile Artery Forceps Straight	14cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Crile Artery Forceps Curved	14cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Halstead-Mosquito Artery Forceps Straight	12.5cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Halstead-Mosquito Artery Forceps Curved	12.5cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Adson Dressing Forcep Standard	12cm	KLINI	1 x	KLINI	1 x	ARMO	1 x	
Adson Brown Tissue Forcep	12cm	KLINI	1 x	KLINI	1 x	ARMO	1 x	
Adson Tissue Forcep 1x2 teeth Standard	12cm	KLINI	1 x	KLINI	1 x	ARMO	1 x	
Allis Tissue Forceps 4x5 teeth	15cm	ARMO	4 x	ARMO	4 x	ARMO	4 x	
Mayo Scissors Straight	14.5cm	HIPP	1 x	KLINI	1 x	ARMO	1 x	
Metzenbaum Scissors Curved blunt/blunt	15.5cm	HIPP	1 x	KLINI	1 x	ARMO	1 x	
Senn Miller Retractor 22x7mm blades Blunt	16cm	ARMO	2 x	ARMO	2 x	ARMO	2 x	
Spencer Scissors	9cm	ARMO	1 x	ARMO	1 x	ARMO	1 x	
Scalpel Handle standard grip	no. 4	ARMO	1 x	ARMO	1 x	ARMO	1 x	
Backhaus Towel Clips	9cm	ARMO	4 x	ARMO	4 x	ARMO	4 x	



Ophthalmic kits



A9021-G1 V A9021-G2 V A9021-G3 V wire

A9021-G3

OPHTHALMIC KITS	A9021-G1 V		A9021-G2 V		A9021-G3 V wire		A9021-G3		
Instrument Inclusions	Brand	Qty	Brand	Qty	Brand	Qty	Brand	Qty	
Adson Dressing Forcep Standard	12cm	KLINI	1 x	KLINI	1 x	ARMO	1 x	ARMO	1 x
Adson Tissue Forceps 1x 2 teeth Micro	12cm	HIPP	1 x	KLINI	1 x	ARMO	1 x	ARMO	1 x
Desmarres Chalazion (No.2)	9cm	ARMO	1 x	ARMO	1 x	ARMO	1 x	ARMO	1 x
Iris Scissor Straight	11cm	HIPP	2 x	KLINI	2 x	ARMO	2 x	ARMO	2 x
Strabismus Scissor Straight	11cm	HIPP	1 x	KLINI	1 x	ARMO	1 x	ARMO	1 x
Castroviejo Ophthalmic Needle Holder (No Lock) Smooth Jaw	13cm	ARMO	1 x	ARMO	1 x	ARMO	1 x	ARMO	1 x
Scalpel Handle No.7 for Blades 10-15	No. 7	ARMO	1 x	ARMO	1 x	ARMO	1 x	ARMO	1 x
Barraquer Wire Adult Eye Specula							VINZENZ SATTLER	1 x	
Lieberman V Wire Eye Specula		VINZENZ SATTLER	1 x	VINZENZ SATTLER	1 x	VINZENZ SATTLER	1 x		



19/9 Hudson Avenue
Castle Hill NSW 2154
Australia
phone (02) 9659 9548
email admin@bydandmedical.com.au

bydandmedical.com.au