


# SURGICRYL<sup>®</sup>

## PGA

### POLYGLYCOLIC ACID

SURGICRYL<sup>®</sup>  
PGA  
POLYGLYCOLIC ACID



<b>Type</b>	Braided and coated multifilament
<b>Composition</b>	Polyglycolic acid
<b>Coating</b>	Polycaprolactone and calcium stearate (1%)
<b>Colour</b>	Violet / Beige - undyed
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. Approximately 50% of tensile strength remain after 21 days.
<b>Presentation</b>	Box with 12 sutures



*For quality  
&  
safety*

**smi**

[www.sutures.be](http://www.sutures.be)

# SURGICRYL<sup>®</sup>

PGA  
POLYGLYCOLIC ACID

## CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Braided structure
  
- Coated multifilament
  
- Hermetically sealed packing

## BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Excellent handling properties
  - High flexibility
  - High tensile strength
  - Excellent knot security
- Soft passage through the tissues with minimal sawing, tissue drag and trauma
- Good knot tie down / adjustability
- Reduced capillary action
  
- Guaranteed seal and product sterility

**INDICATIONS** General soft tissue closing and/or ligation especially in


- General surgery
- Skin closure
- Gastrointestinal surgery
- Gynaecology
- Obstetrics
- Plastic surgery
- Urology
- Ophthalmic surgery
- Orthopaedics

**CONTRAINDICATIONS** Cardiovascular surgery  
Neurosurgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

## TENSILE STRENGTH



## ABSORPTION



Without needle



Metric	USP	Length		Order number
2	3/0	150 cm	<i>violet</i>	11200015
3	2/0	150 cm	<i>violet</i>	11300015
3.5	0	150 cm	<i>violet</i>	11350015
4	1	150 cm	<i>violet</i>	11400015
5	2	150 cm	<i>violet</i>	11500015
6	3&4	150 cm	<i>violet</i>	11600015
7	5	150 cm	<i>violet</i>	11700015
8	6	150 cm	<i>violet</i>	11800015
1.5	4/0	3 x 50 cm	<i>violet</i>	11150050
2	3/0	3 x 50 cm	<i>violet</i>	11200050
3	2/0	3 x 50 cm	<i>violet</i>	11300050
3.5	0	3 x 50 cm	<i>violet</i>	11350050
4	1	3 x 50 cm	<i>violet</i>	11400050
5	2	3 x 50 cm	<i>violet</i>	11500050
2	3/0	12 x 45 cm	<i>violet</i>	11200045
3	2/0	12 x 45 cm	<i>violet</i>	11300045
3.5	0	12 x 45 cm	<i>violet</i>	11350045
4	1	12 x 45 cm	<i>violet</i>	11400045
5	2	12 x 45 cm	<i>violet</i>	11500045
1.5	4/0	6 x 45 cm	<i>violet</i>	11150045-6
2	3/0	6 x 45 cm	<i>violet</i>	11200045-6
3	2/0	6 x 45 cm	<i>violet</i>	11300045-6
3.5	0	6 x 45 cm	<i>violet</i>	11350045-6
4	1	6 x 45 cm	<i>violet</i>	11400045-6
5	2	6 x 45 cm	<i>violet</i>	11500045-6
2	3/0	10 x 60 cm	<i>violet</i>	11200060
3	2/0	10 x 60 cm	<i>violet</i>	11300060
3.5	0	10 x 60 cm	<i>violet</i>	11350060
2	3/0	5 x 70 cm	<i>violet</i>	11200070
3	2/0	5 x 70 cm	<i>violet</i>	11300070
3.5	0	5 x 70 cm	<i>violet</i>	11350070
4	1	5 x 70 cm	<i>violet</i>	11400070
5	2	5 x 70 cm	<i>violet</i>	11500070

## Without needle



### Ligapack



Metric	USP	Length	Order number
1.5	4/0	250 cm	<i>violet</i> 11150025
2	3/0	250 cm	<i>violet</i> 11200025
3	2/0	250 cm	<i>violet</i> 11300025
3.5	0	250 cm	<i>violet</i> 11350025
4	1	250 cm	<i>violet</i> 11400025
5	2	250 cm	<i>violet</i> 11500025
6	3&4	250 cm	<i>violet</i> 11600025














## Round bodied taper point










### 1/2 Circle













	Metric	USP	Length	Order number
HR 13 mm	0.7	6/0	75 cm	<i>violet</i> 11070113
	1	5/0	75 cm	<i>violet</i> 11100113
	1.5	4/0	75 cm	<i>violet</i> 11150113
HR 17 mm	0.7	6/0	75 cm	<i>violet</i> 11070117
	1	5/0	75 cm	<i>violet</i> 11100117
	1.5	4/0	75 cm	<i>violet</i> 11150117
	2	3/0	75 cm	<i>violet</i> 11200117
	3	2/0	75 cm	<i>violet</i> 11300117
HR 22 mm	1.5	4/0	75 cm	<i>violet</i> 11150122
	2	3/0	75 cm	<i>violet</i> 11200122
	2	3/0	75 cm	<i>undyed</i> 11200122U
HR 26 mm	3	2/0	75 cm	<i>violet</i> 11300122
	1.5	4/0	75 cm	<i>violet</i> 11150126
	2	3/0	75 cm	<i>violet</i> 11200126
HR 30 mm	3	2/0	75 cm	<i>violet</i> 11300126
	3	2/0	75 cm	<i>undyed</i> 11300126U
	3.5	0	75 cm	<i>violet</i> 11350126
	4	1	75 cm	<i>violet</i> 11400126
	2	3/0	75 cm	<i>violet</i> 11200130
HR 36 mm	3	2/0	75 cm	<i>violet</i> 11300130
	3	2/0	75 cm	<i>undyed</i> 11300130U
	3.5	0	75 cm	<i>violet</i> 11350130
	4	1	75 cm	<i>violet</i> 11400130
	2	3/0	75 cm	<i>violet</i> 11200136
HR 36 mm	3	2/0	75 cm	<i>violet</i> 11300136
	3.5	0	90 cm	<i>violet</i> 11350136
	4	1	90 cm	<i>violet</i> 11400136

Round bodied taper point			Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HR 40 mm		3	2/0	75 cm	<i>violet</i> 11300140
			3.5	0	90 cm	<i>violet</i> 11350140
			4	1	75 cm	<i>violet</i> 11400140-75
			4	1	90 cm	<i>violet</i> 11400140
			5	2	75 cm	<i>violet</i> 11500140-75
	HR 48 mm		5	2	90 cm	<i>violet</i> 11500140
			3	2/0	75 cm	<i>violet</i> 11300148
			3.5	0	90 cm	<i>violet</i> 11350148
			4	1	90 cm	<i>violet</i> 11400148
			5	2	90 cm	<i>violet</i> 11500148
6	3&4	90 cm	<i>violet</i> 11600148			
Round bodied taper point			Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DR 17 mm		1.5	4/0	75 cm	<i>violet</i> 11150217
			2	3/0	75 cm	<i>violet</i> 11200217
Round bodied taper point			Metric	USP	Length	Order number
<b>5/8 Circle</b> 	FR 27 mm		3	2/0	75 cm	<i>violet</i> 11300327
Round bodied taper point			Metric	USP	Length	Order number
<b>J-needle</b> 	JR 36 mm		3	2/0	75 cm	<i>violet</i> 11303036
			3.5	0	75 cm	<i>violet</i> 11353036
			4	1	75 cm	<i>violet</i> 11403040

Round bodied taper point		Metric	USP	Length	Order number
5/8 Circle  J-needle 	JR 36 mm FR 27 mm 	3.5	0	75 cm	<i>violet</i> 11353036-0327

Round bodied taper point loop		Metric	USP	Length	Order number
1/2 Circle 	HR 50 mm 	4	1	150 cm	<i>violet</i> 11400950

Round bodied taper cutting point		Metric	USP	Length	Order number
1/2 Circle 	HRT 26 mm 	3	2/0	75 cm	<i>violet</i> 11300526
	HRT 30 mm 	3	2/0	75 cm	<i>violet</i> 11300530
		3.5	0	75 cm	<i>violet</i> 11350530
		4	1	75 cm	<i>violet</i> 11400530
	HRT 36 mm 	3	2/0	90 cm	<i>violet</i> 11300536
		3.5	0	90 cm	<i>violet</i> 11350536
		4	1	90 cm	<i>violet</i> 11400536
	HRT 40 mm 	3.5	0	90 cm	<i>violet</i> 11350540
		4	1	90 cm	<i>violet</i> 11400540
	HRT 48 mm 	5	2	90 cm	<i>violet</i> 11500540
HRT 48 mm	4	1	90 cm	<i>violet</i> 11400548	

Round bodied taper cutting point		Metric	USP	Length	Order number
J-needle 	JRT 31 mm 	4	1	75 cm	<i>violet</i> 11403131

Reverse cutting













3/8 Circle





		Metric	USP	Length		Order number
DS 12 mm		0.7	6/0	45 cm	<i>violet</i>	11071512
		0.7	6/0	45 cm	<i>undyed</i>	11071512U
		1	5/0	45 cm	<i>violet</i>	11101512
		1.5	4/0	75 cm	<i>violet</i>	11151512
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>undyed</i>	11071512UP
		1	5/0	50 cm	<i>undyed</i>	11101512UP
		1.5	4/0	50 cm	<i>undyed</i>	11151512UP
DS 16 mm		0.7	6/0	75 cm	<i>violet</i>	11071516
		1	5/0	75 cm	<i>violet</i>	11101516
		1	5/0	75 cm	<i>undyed</i>	11101516U
		1.5	4/0	75 cm	<i>violet</i>	11151516
		1.5	4/0	75 cm	<i>undyed</i>	11151516U
		2	3/0	75 cm	<i>violet</i>	11201516
DS 16 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151516UP
DS 19 mm		1	5/0	75 cm	<i>violet</i>	11101519
		1.5	4/0	75 cm	<i>violet</i>	11151519
		1.5	4/0	75 cm	<i>undyed</i>	11151519U
		2	3/0	75 cm	<i>violet</i>	11201519
		2	3/0	75 cm	<i>undyed</i>	11201519U
		3	2/0	75 cm	<i>violet</i>	11301519
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i>	11151519UP
		2	3/0	50 cm	<i>undyed</i>	11201519UP
DS 24 mm		2	3/0	75 cm	<i>violet</i>	11201524
		2	3/0	75 cm	<i>undyed</i>	11201524U
		3	2/0	75 cm	<i>violet</i>	11301524
		3	2/0	75 cm	<i>undyed</i>	11301524U
		3.5	0	75 cm	<i>violet</i>	11351524
DS 24 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	11201524UP
		3	2/0	50 cm	<i>undyed</i>	11301524UP
DS 30 mm		1.5	4/0	75 cm	<i>violet</i>	11151530
		2	3/0	75 cm	<i>violet</i>	11201530
		3	2/0	75 cm	<i>violet</i>	11301530
		3	2/0	75 cm	<i>undyed</i>	11301530U
		3.5	0	75 cm	<i>violet</i>	11351530
		4	1	75 cm	<i>violet</i>	11401530
		5	2	75 cm	<i>violet</i>	11501530







Reverse cutting		Metric	USP	Length	Order number
3/8 Circle 	DS 36 mm 	3	2/0	75 cm	<i>violet</i> 11301536
		3.5	0	75 cm	<i>violet</i> 11351536
		4	1	75 cm	<i>violet</i> 11401536
		5	2	75 cm	<i>violet</i> 11501536
		3	2/0	75 cm	<i>violet</i> 11301540
	DS 40 mm 	3.5	0	75 cm	<i>violet</i> 11351540
		4	1	75 cm	<i>violet</i> 11401540
		5	2	75 cm	<i>violet</i> 11501540



Reverse cutting		Metric	USP	Length	Order number
1/2 Circle 	HS 22 mm 	2	3/0	75 cm	<i>violet</i> 11201422
		2	3/0	75 cm	<i>undyed</i> 11201422U
		3	2/0	75 cm	<i>violet</i> 11301422
		3	2/0	75 cm	<i>undyed</i> 11301422U
	HS 26 mm 	1.5	4/0	75 cm	<i>violet</i> 11151426
		2	3/0	75 cm	<i>violet</i> 11201426
		3	2/0	75 cm	<i>violet</i> 11301426
		3.5	0	75 cm	<i>violet</i> 11351426
	HS 30 mm 	5	2	90 cm	<i>violet</i> 11501430
		3	2/0	90 cm	<i>violet</i> 11301436
		3.5	0	90 cm	<i>violet</i> 11351436
		4	1	90 cm	<i>violet</i> 11401436
	HS 36 mm 	5	2	90 cm	<i>violet</i> 11501436
		3	2/0	90 cm	<i>violet</i> 11301440
		3.5	0	90 cm	<i>violet</i> 11351440
		4	1	90 cm	<i>violet</i> 11401440
	HS 40 mm 	5	2	90 cm	<i>violet</i> 11501440
		3.5	0	90 cm	<i>violet</i> 11351448
		4	1	90 cm	<i>violet</i> 11401448
	HS 48 mm 	5	2	90 cm	<i>violet</i> 11501448
6		3&4	90 cm	<i>violet</i> 11601448	







Cutting		Metric	USP	Length		Order number	
Straight		GS 60 mm	2	3/0	75 cm	<i>violet</i>	11202160
			2	3/0	75 cm	<i>undyed</i>	11202160U
			3	2/0	75 cm	<i>violet</i>	11302160

Cutting		Metric	USP	Length		Order number	
3/8 Circle		NS 16 mm	1	5/0	45 cm	<i>violet</i>	11101716

Side cutting spatula double armed		Metric	USP	Length		Order number	
1/2 Circle		2 x HSP 8 mm 200 µm	0.5	7/0	45 cm	<i>violet</i>	OSU07

Side cutting spatula double armed		Metric	USP	Length		Order number	
3/8 Circle		2 x DSP 6 mm 200 µm	0.4	8/0	30 cm	<i>violet</i>	OSU08-30
			0.4	8/0	45 cm	<i>violet</i>	OSU08
		2 x DSP 6,2 mm 150 µm	0.2	10/0	30 cm	<i>violet</i>	OSU10
		2 x DSP 6,4 mm 200 µm	0.3	9/0	30 cm	<i>violet</i>	OSU09
			0.5	7/0	45 cm	<i>violet</i>	OSU072
		2 x DSP 8 mm 200 µm	0.5	7/0	45 cm	<i>violet</i>	OSU071

Side cutting spatula		Metric	USP	Length		Order number	
3/8 circle		DSP 6,2 mm 150 µm	0.2	10/0	10 cm	<i>violet</i>	OSU101
		DSP 6,4 mm 200 µm	0.3	9/0	30 cm	<i>violet</i>	OSU09-1

Side cutting spatula double armed		Metric	USP	Length		Order number	
1/4 Circle		2 x VSP 8 mm 300 µm 350 µm	0.7	6/0	45 cm	<i>violet</i>	OSU06
			1	5/0	45 cm	<i>violet</i>	OSU05

## IN CASSETTES

Use within 4 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	11151
2	3/0	15 m	11201
2	3/0	25 m	11202
3	2/0	15 m	11301
3	2/0	25 m	11302
3	2/0	50 m	11305
3.5	0	15 m	11351
3.5	0	25 m	11352
3.5	0	50 m	11355
4	1	15 m	11401
4	1	25 m	11402
4	1	50 m	11405
5	2	15 m	11501
5	2	25 m	11502
5	2	50 m	11505
5	2	100 m	115010
6	3&4	15 m	11601
6	3&4	25 m	11602
6	3&4	50 m	11605
6	3&4	100 m	116010
7	5	15 m	11701
7	5	25 m	11702
7	5	50 m	11705
8	6	15 m	11801
8	6	25 m	11802
8	6	50 m	11805



## with storage solution

Metric	USP	Length	Order number
1.5	4/0	15 m	15151
2	3/0	15 m	15201
2	3/0	25 m	15202
3	2/0	15 m	15301
3	2/0	25 m	15302
3	2/0	50 m	15305
3.5	0	15 m	15351
3.5	0	25 m	15352
3.5	0	50 m	15355
4	1	15 m	15401
4	1	25 m	15402
4	1	50 m	15405
5	2	15 m	15501

Metric	USP	Length	Order number
5	2	25 m	15502
5	2	50 m	15505
5	2	100 m	155010
6	3&4	15 m	15601
6	3&4	25 m	15602
6	3&4	50 m	15605
6	3&4	100 m	156010
7	5	15 m	15701
7	5	25 m	15702
7	5	50 m	15705
8	6	15 m	15801
8	6	25 m	15802
8	6	50 m	15805

Use within 6 months after first opening