


# SURGICRYL<sup>®</sup>

## MONOFAST

### POLYGLCAPRONE 25



<b>Type</b>	Monofilament
<b>Composition</b>	Polyglcaprone 25, a copolymer made of 75% glycolide and 25% $\epsilon$ -caprolactone
<b>Coating</b>	None
<b>Colour</b>	Violet / Beige - undyed
<b>Tissue reaction</b>	Minimal
<b>Absorption</b>	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. Approximately 60% of tensile strength remain after 7 days.
<b>Presentation</b>	Box with 12 sutures



*For quality  
&  
safety*

The logo for smi, consisting of the lowercase letters 'smi' in a bold, blue, sans-serif font with a white outline.

[www.sutures.be](http://www.sutures.be)

# SURGICRYL<sup>®</sup>

## MONOFAST

### POLYGLECAPRONE 25

#### CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter
  
- Pliability and softness
- Hermetically sealed packing

#### BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
  
- Good handling properties
- Easy and secure knotting
  
- Guaranteed seal and product sterility


**INDICATIONS** Superficial soft tissue and/or ligation  
 especially in General surgery  
 Skin closure  
 Plastic surgery  
 Gastrointestinal surgery  
 Gynaecology  
 Urology

**CONTRAINDICATIONS** Cardiovascular surgery  
 Neurosurgery  
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

#### TENSILE STRENGTH



#### ABSORPTION



Without needle 

**Metric**    **USP**    **Length**    **Order number**

3	2/0	150 cm	<i>violet</i>	16300015
3.5	0	150 cm	<i>violet</i>	16350015
4	1	150 cm	<i>violet</i>	16400015
5	2	120 cm	<i>violet</i>	16500012
5	2	150 cm	<i>violet</i>	16500015

Without needle 

**Metric**    **USP**    **Length**    **Order number**







**Ligapack**



5	2	250 cm	<i>violet</i>	16500025
---	---	--------	---------------	----------

Round bodied taper point 

**Metric**    **USP**    **Length**    **Order number**

<b>1/2 Circle</b> 	HR 13 mm		0.7	6/0	45 cm	<i>violet</i>	16070113
			1	5/0	45 cm	<i>violet</i>	16100113
			1.5	4/0	75 cm	<i>violet</i>	16150113
	HR 17 mm		1	5/0	45 cm	<i>violet</i>	16100117
			1.5	4/0	75 cm	<i>violet</i>	16150117
			2	3/0	75 cm	<i>violet</i>	16200117
	HR 22 mm		1.5	4/0	75 cm	<i>violet</i>	16150122
			2	3/0	75 cm	<i>violet</i>	16200122
	HR 26 mm		1.5	4/0	75 cm	<i>violet</i>	16150126
			2	3/0	75 cm	<i>violet</i>	16200126
			3	2/0	75 cm	<i>violet</i>	16300126
			3.5	0	75 cm	<i>violet</i>	16350126
HR 30 mm		4	1	75 cm	<i>violet</i>	16400126	
		2	3/0	75 cm	<i>violet</i>	16200130	
		3	2/0	75 cm	<i>violet</i>	16300130	
		3.5	0	75 cm	<i>violet</i>	16350130	
			4	1	75 cm	<i>violet</i>	16400130

**Round bodied taper point**



**1/2 Circle**



		Metric	USP	Length		Order number	
HR 36 mm		3	2/0	75 cm	<i>violet</i>	16300136	
		3.5	0	90 cm	<i>violet</i>	16350136	
		4	1	90 cm	<i>violet</i>	16400136	
	HR 40 mm		3	2/0	75 cm	<i>violet</i>	16300140
			3.5	0	90 cm	<i>violet</i>	16350140
			4	1	90 cm	<i>violet</i>	16400140
HR 48 mm		3.5	0	90 cm	<i>violet</i>	16350148	
		4	1	90 cm	<i>violet</i>	16400148	





**Reverse cutting**










**3/8 Circle**







		Metric	USP	Length		Order number
DS 12 mm		0.7	6/0	45 cm	<i>undyed</i>	16071512
		1	5/0	45 cm	<i>undyed</i>	16101512
		1	5/0	45 cm	<i>violet</i>	16101512V
DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>undyed</i>	16071512P
		1	5/0	50 cm	<i>undyed</i>	16101512P
		1.5	4/0	50 cm	<i>undyed</i>	16151512P
DS 16 mm		1	5/0	75 cm	<i>undyed</i>	16101516
		1.5	4/0	75 cm	<i>undyed</i>	16151516
		2	3/0	75 cm	<i>undyed</i>	16201516
DS 16 mm Easy Pass		1	5/0	50 cm	<i>undyed</i>	16101516P
		1.5	4/0	50 cm	<i>undyed</i>	16151516P
		1	5/0	45 cm	<i>undyed</i>	16101519
DS 19 mm		1.5	4/0	75 cm	<i>undyed</i>	16151519
		1.5	4/0	75 cm	<i>violet</i>	16151519V
		2	3/0	75 cm	<i>undyed</i>	16201519
		3	2/0	75 cm	<i>undyed</i>	16301519
		1.5	4/0	50 cm	<i>undyed</i>	16151519P
DS 19 mm Easy Pass		2	3/0	50 cm	<i>undyed</i>	16201519P
		2	3/0	75 cm	<i>undyed</i>	16201524
DS 24 mm		3	2/0	75 cm	<i>undyed</i>	16301524

Reverse cutting 		Metric	USP	Length	Order number
<b>3/8 Circle</b> 	DS 26 mm 	1.5	4/0	75 cm	<i>undyed</i> 16151526
		2	3/0	75 cm	<i>undyed</i> 16201526
		2	3/0	75 cm	<i>violet</i> 16201526V
	DS 40 mm 	3	2/0	75 cm	<i>undyed</i> 16301526
		3	2/0	75 cm	<i>violet</i> 16301526V
		3.5	0	75 cm	<i>undyed</i> 16351540
		4	1	75 cm	<i>undyed</i> 16401540

Reverse cutting 		Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HS 36 mm 	3	2/0	75 cm	<i>undyed</i> 16301436
		3.5	0	75 cm	<i>undyed</i> 16351436
	HS 40 mm 	4	1	90 cm	<i>undyed</i> 16401440

Cutting 		Metric	USP	Length	Order number
<b>Straight</b> 	GS 60 mm 	2	3/0	75 cm	<i>undyed</i> 16202160
		3	2/0	75 cm	<i>undyed</i> 16302160

Blunt 		Metric	USP	Length	Order number
<b>1/2 Circle</b> 	HRB 40 mm 	3.5	0	90 cm	<i>violet</i> 16351840
	HRB 48 mm 	4	1	90 cm	<i>violet</i> 16401848