


SILK



	<i>Silk Braided</i>	<i>Silk Virgin</i>
Type	Braided multifilament	Twisted multifilament
Composition	Braided fibres from the cocoon of the silkworm	Twisted fibres from the cocoon of the silkworm
Coating	Waxes Silicone	None
Colour	Black	Blue
Tissue reaction	Moderate	
Absorption	Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.	
Presentation	Box with 12 sutures	

*For quality
&
safety*

smi

www.sutures.be

SILK

CHARACTERISTICS

- ☑ High tensile strength, non absorbable
- ☑ Braided or twisted structure

- ☑ Coated multifilament

- ☑ Hermetically sealed packing

BENEFITS

- Good and extended tissue support for up to 3 month

- Excellent handling properties
High flexibility
High tensile strength
Excellent knot security

- Soft passage through the tissues with minimal sawing,
tissue drag and trauma
Good knot tie down / adjustability
Reduced capillary action

- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing
and/or ligation
especially in General surgery
Skin closure
Ophthalmic surgery
Oral surgery
Gastrointestinal surgery
Gynaecology
Obstetrics

CONTRAINDICATIONS Urinary tract
Biliary tract

Known allergies or sensitivities to silk

Without needle



Metric **USP** **Length** **Order number**

2	3/0	150 cm	<i>black</i>	8200015
3	2/0	150 cm	<i>black</i>	8300015
3.5	0	150 cm	<i>black</i>	8350015
4	1	150 cm	<i>black</i>	8400015
5	2	150 cm	<i>black</i>	8500015
6	3&4	150 cm	<i>black</i>	860 0015
8	6	150 cm	<i>black</i>	8800015
2	3/0	6 x 45 cm	<i>black</i>	8200044
3	2/0	6 x 45 cm	<i>black</i>	8300044
3.5	0	6 x 45 cm	<i>black</i>	8350044
4	1	6 x 45 cm	<i>black</i>	8400044
5	2	6 x 45 cm	<i>black</i>	8500044
1.5	4/0	17 x 45 cm	<i>black</i>	8150046
2	3/0	17 x 45 cm	<i>black</i>	8200046
3	2/0	17 x 45 cm	<i>black</i>	8300046
3.5	0	17 x 45 cm	<i>black</i>	8350046
4	1	17 x 45 cm	<i>black</i>	8400046
5	2	17 x 45 cm	<i>black</i>	8500046
1.5	4/0	13 x 60 cm	<i>black</i>	8150060
2	3/0	13 x 60 cm	<i>black</i>	8200060
3	2/0	13 x 60 cm	<i>black</i>	8300060
3.5	0	13 x 60 cm	<i>black</i>	8350060
4	1	13 x 60 cm	<i>black</i>	8400060
2	3/0	12 x 75 cm	<i>black</i>	8200075
3	2/0	12 x 75 cm	<i>black</i>	8300075
3.5	0	12 x 75 cm	<i>black</i>	8350075
4	1	12 x 75 cm	<i>black</i>	8400075
5	2	12 x 75 cm	<i>black</i>	8500075

Without needle



Ligapack












Metric **USP** **Length** **Order number**

2	3/0	250 cm	<i>black</i>	8200025
3	2/0	250 cm	<i>black</i>	8300025
3.5	0	250 cm	<i>black</i>	8350025
4	1	250 cm	<i>black</i>	8400025
5	2	250 cm	<i>black</i>	8500025
6	3&4	250 cm	<i>black</i>	860 0025






Round bodied taper point

1/2 Circle




		Metric	USP	Length	Color	Order number
HR 12 mm		0.7	6/0	45 cm	black	8070112
		1	5/0	45 cm	black	8100112
		1.5	4/0	75 cm	black	8150112
HR 16 mm		0.7	6/0	45 cm	black	8070116
		1	5/0	75 cm	black	8100116
		1.5	4/0	75 cm	black	8150116
HR 20 mm		2	3/0	75 cm	black	8200116
		1	5/0	75 cm	black	8100120
		1.5	4/0	75 cm	black	8150120
HR 26 mm		2	3/0	75 cm	black	8200120
		3	2/0	75 cm	black	8300120
		3.5	0	75 cm	black	8350120
		1	5/0	75 cm	black	8100126
		1.5	4/0	75 cm	black	8150126
HR 30 mm		2	3/0	75 cm	black	8200126
		3	2/0	75 cm	black	8300126
		3.5	0	75 cm	black	8350126
		4	1	75 cm	black	8400126
		5	2	75 cm	black	8500126
		2	3/0	75 cm	black	8200130
		3	2/0	75 cm	black	8300130
HR 37 mm		3.5	0	75 cm	black	8350130
		4	1	75 cm	black	8400130
		5	2	75 cm	black	8500130
		2	3/0	75 cm	black	8200137
HR 40 mm		3	2/0	75 cm	black	8300137
		3.5	0	75 cm	black	8350137
		4	1	75 cm	black	8400137
		5	2	75 cm	black	8500137
HR 45 mm		4	1	75 cm	black	8400140
		5	2	75 cm	black	8500140
HR 45 mm		3.5	0	75 cm	black	8350145
		4	1	75 cm	black	8400145
		5	2	75 cm	black	8500145







Round bodied taper point

			Metric	USP	Length	Color	Order number
3/8 Circle 	DR 12 mm 	1	5/0	45 cm	<i>black</i>	8100212	
		1.5	4/0	75 cm	<i>black</i>	8150216	
	DR 16 mm 	2	3/0	75 cm	<i>black</i>	8200216	
		2	3/0	75 cm	<i>black</i>	8200225	
		3	2/0	75 cm	<i>black</i>	8300225	
		3.5	0	75 cm	<i>black</i>	8350225	
	DR 25 mm 	4	1	75 cm	<i>black</i>	8400225	
		3	2/0	75 cm	<i>black</i>	8300230	
	DR 30 mm 						

Round bodied taper point










		Metric	USP	Length	Color	Order number
Straight	GR 40 mm 	3	2/0	75 cm	<i>black</i>	8302040

Reverse cutting

			Metric	USP	Length	Color	Order number
3/8 Circle 	DS 12 mm 	0.5	7/0	45 cm	<i>black</i>	8051512	
		0.7	6/0	45 cm	<i>black</i>	8071512	
		1	5/0	45 cm	<i>black</i>	8101512	
		1.5	4/0	75 cm	<i>black</i>	8151512	
		2	3/0	75 cm	<i>black</i>	8201512	
	DS 12 mm Easy Pass 	0.7	6/0	50 cm	<i>black</i>	8071512P	
	DS 14 mm 	0.7	6/0	45 cm	<i>black</i>	8071514	
	DS 16 mm 	0.7	6/0	75 cm	<i>black</i>	8071516	
		1	5/0	75 cm	<i>black</i>	8101516	
		1.5	4/0	75 cm	<i>black</i>	8151516	
		2	3/0	75 cm	<i>black</i>	8201516	
		3	2/0	75 cm	<i>black</i>	8301516	
	DS 16 mm Easy Pass 	1	5/0	50 cm	<i>black</i>	8101516P	
		1.5	4/0	50 cm	<i>black</i>	8151516P	


Reverse cutting



		Metric	USP	Length	Color	Order number
DS 19 mm		0.7	6/0	75 cm	black	8071519
		1	5/0	75 cm	black	8101519
		1.5	4/0	75 cm	black	8151519
		2	3/0	75 cm	black	8201519
		3	2/0	75 cm	black	8301519
		3.5	0	75 cm	black	8351519
		4	1	75 cm	black	8401519
DS 19 mm Easy Pass		1.5	4/0	50 cm	black	8151519P
		2	3/0	50 cm	black	8201519P
DS 24 mm		1	5/0	75 cm	black	8101524
		1.5	4/0	75 cm	black	8151524
		2	3/0	75 cm	black	8201524
		3	2/0	75 cm	black	8301524
		3.5	0	75 cm	black	8351524
		4	1	75 cm	black	8401524
DS 24 mm Easy Pass		2	3/0	75 cm	black	8201524P
DS 26 mm		2	3/0	75 cm	black	8201526
		3	2/0	75 cm	black	8301526
DS 30 mm		1	5/0	75 cm	black	8101530
		1.5	4/0	75 cm	black	8151530
		2	3/0	75 cm	black	8201530
		3	2/0	75 cm	black	8301530
		3.5	0	75 cm	black	8351530
		4	1	75 cm	black	8401530
		5	2	75 cm	black	8501530
DS 37 mm		2	3/0	75 cm	black	8201537
		3	2/0	75 cm	black	8301537
		3.5	0	75 cm	black	8351537
		4	1	75 cm	black	8401537
		5	2	75 cm	black	8501537
		6	3&4	75 cm	black	8601537
DS 45 m		3	2/0	75 cm	black	8301545
		3.5	0	75 cm	black	8351545
		4	1	75 cm	black	8401545
		5	2	75 cm	black	8501545
DS 60 mm		3	2/0	75 cm	black	8301560
		4	1	75 cm	black	8401560

3/8 Circle



Reverse cutting  **Metric** **USP** **Length** **Order number**

3/8 Circle



DS 90 mm



3.5

0

100 cm

black

8351590

Reverse cutting  **Metric** **USP** **Length** **Order number**

1/2 Circle



HS 12 mm



0.7

6/0

75 cm

black

8071412

HS 17 mm



1

5/0

75 cm

black

8101417

1.5

4/0

75 cm

black

8151417

2

3/0

75 cm

black

8201417

HS 20 mm



1.5

4/0

75 cm

black

8151420

2

3/0

75 cm

black

8201420

HS 22 mm



1

5/0

75 cm

black

8101422

1.5

4/0

75 cm

black

8151422

2

3/0

75 cm

black

8201422

3

2/0

75 cm

black

8301422

HS 26 mm



1

5/0

75 cm

black

8101426

1.5

4/0

75 cm

black

8151426

2

3/0

75 cm

black

8201426

3

2/0

75 cm

black

8301426

3.5

0

75 cm

black

8351426

HS 36 mm



3.5

0

75 cm

black

8351436

4

1

75 cm

black

8401436

Cutting  **Metric** **USP** **Length** **Order number**

Straight

GS 40 mm



3

2/0

75 cm

black

8302140

2

3/0

75 cm

black

8202160

GS 60 mm



3

2/0

75 cm

black

8302160

3.5

0

75 cm

black

835 2160

4

1

75 cm

black

840 2160





5




2

75 cm

black

850 2160

Round bodied tapercut point 		Metric	USP	Length	Order number
1/2 Circle 	HRT 18 mm 	1.5	4/0	75 cm	<i>black</i> 8150518
		2	3/0	75 cm	<i>black</i> 8200518
	HRT 27 mm 	2	3/0	75 cm	<i>black</i> 8200527
		3	2/0	75 cm	<i>black</i> 8300527

Round bodied tapercut point 		Metric	USP	Length	Order number
3/8 Circle 	DRT 18 mm 	1.5	4/0	75 cm	<i>black</i> 8150418
		2	3/0	75 cm	<i>black</i> 8200418